

# biocap

Needle-free drug aspiration; for multiple vial access

- Contains Latex: No
- Contains DEHP: No
- Contains biological or animal-based product: No
- Non-pyrogenic: Yes

Product Code	Cap Ø	Units/Box
821.21	14 mm	50
821.22	20 mm	50



# bionector

Bionector is a multi-purpose, closed, needle-free connection device.

When Bionector is in place, you can infuse, inject sample and change the I.V. tubing without the I.V. circuit to the atmosphere. Bionector's membrane automatically opens the fluid pathway only when a male Luer is disconnected, the membrane automatically seals the fluid pathway.

- Bionector preserves constantly a closed system
- 7 days or 360 connections
- Straight fluid pathway
- Neutral displacement of fluid
- Minimal dead space
- Can be used with power injectors (CT-scan).  
Pressure resistance: 350 psi
- MRI conditional
- Contains Latex: No
- Contains DEHP: No
- Contains biological or animal-based product: No
- Non-pyrogenic: Yes



Product Code	Flow rate	Priming volume	Equivalent Gauge G	Number of accesses	Max. duration of use	Length	Diameter	Units/Box
896.038	105 ml/min	0.03 ml	-	360	7 days	25.7 mm	10 mm	50