

COMFIT™ Endotracheal Tube Holder

The Tapeless Way to Secure an Endotracheal Tube



COMFIT™ is a reliable, tapeless solution used to secure endotracheal tubes and better manage patient airways. Physicians appreciate COMFIT'S compliance with Advanced Cardiac Life Support (ACLS) Guidelines that recommend endotracheal tube stabilization devices as a substitute for tape. For patients, COMFIT™ provides a comfortable fit without the skin irritation that tape can cause.

Features and Benefits

- Significantly decreases the incidence of self-extubation¹
- Wide cotton-lined neckband minimizes skin irritation, providing maximum patient comfort
- No tape on the tube or the patient
- Completely adjustable
- Minimal plastic loop around the endotracheal tube allows access to the oral cavity
- A low-cost device that requires infrequent changing
- Latex-free product
- Complies with American Heart Association's ACLS Guidelines recommending the use of commercially available endotracheal tube holders instead of tape
- Plastic double securing loop, ½- x 19-inch neckband



¹ McAndrew J, RR. Reduction of self extubation using Comfit™ tracheal tube holder. Respiratory Care. November 1994; 39 (11).

Item #	Description
42-1820	Tapeless Endotracheal Tube Holder (25/Box) Each set includes: (1) plastic tube grip, with hook/loop strap, foam foot pads, (2) foam "lollipop" tapes

