

dual-flowgastric

Feeding Tubes

The dual-flow gastric tubes are used to prevent the accumulation of gastric secretions and air; thus allowing decompression of the stomach by continuous aspiration without trauma to the gastric mucosa.



Dual-flow tube

Dual-flow gastric tubes

Product Code		Size	Length	Number of eyes	Units/Case
PVC	PUR				
340.06	I340.06	6 Fr	60 cm	6	10
340.08	I340.08	8 Fr	60 cm	6	10
340.10	I340.10	10 Fr	90 cm	10	10
340.12	I340.12	12 Fr	90 cm	10	10
340.14	-	14 Fr	120 cm	10	10
340.16	I340.16	16 Fr	120 cm	10	10
340.18	I340.18	18 Fr	120 cm	10	10