


endotracheal tubes

Single and double lumen tubes in PVC (DEHP-free) for emergency oral or nasal intubation and routine anesthesia. Double lumen tubes have secondary lumen for surfactant administration or airway monitoring.

- Radiopaque
- 15mm adapter
- Markings every half centimeter

Single lumen tubes are also available as a soft tube designed to minimize tracheal damage during long term intubation.



Plain endotracheal tubes

Product Code		Ext. ø	Int. ø	Length	Units/Case
Standard	Soft tube				
516.20	520.20	3.4 mm	2.0 mm	165 mm	25
516.25	520.25	4.1 mm	2.5 mm	165 mm	25
516.30	520.30	4.6 mm	3.0 mm	165 mm	25
516.35	520.35	5.2 mm	3.5 mm	165 mm	25
516.40	520.40	5.7 mm	4.0 mm	230 mm	20
516.45	520.45	6.2 mm	4.5 mm	230 mm	20

Endotracheal tubes with secondary lumen

Product Code	Tube			Secondary lumen		Units/Case
	Ext. ø	Int. ø	Length	Lumen diam.	Flow rate	
5516.20	3.4 mm	2.0 mm	165 mm	0.40 mm	2.15 ml/min	20
5516.25	4.1 mm	2.5 mm	165 mm	0.55 mm	3.35 ml/min	20
5516.30	4.6 mm	3.0 mm	165 mm	0.65 mm	5.00 ml/min	20
5516.35	5.2 mm	3.5 mm	165 mm	0.80 mm	15.0 ml/min	20
5516.40	5.7 mm	4.0 mm	230 mm	0.95 mm	35.0 ml/min	20