

# nutrifit


The perfect fit for pediatrics and adults

Enteral nutrition is recommended when oral food intake proves impossible or insufficient to meet the patient's nutritional needs. Enteral nutrition contributes to:

- ✓ Preserving normal physiology<sup>1</sup>, digestion, and nutrient absorption
- ↘ Reducing:
  - Length of stay<sup>2</sup>
  - Infectious complications<sup>3</sup>
  - Costs<sup>4</sup>
- ⊕ The reinforcement of the immune system<sup>5</sup>



nutrifit  
The ENFit™ brand of Vygon



## PVC FEEDING TUBES

Adapted for children and adults, the PVC nutritif feeding tubes are intended to be used for short-term enteral feeding (up to 5 days) and feature the following characteristics:

- Made of PVC
- Closed distal tip with 2 lateral eyes (4Fr to 12Fr) and 4 lateral eyes (14Fr to 21Fr)
- Small cap with color code to quickly identify the size of the tube
- Numerical markings every centimeter:
  - from 5cm to 35cm for codes 364, 361 & 363
  - from 5cm to 70cm for codes 366, 362
- Radiopaque
- Transparent to better visualize the administered solution.

Code	PVC tube					Color code	Quantity
	Length cm	Fr	Ext. Ø mm	Int. Ø mm	Prim. vol. ml		Units/case
364.047	40	4	1.5	0.8	0.27	Purple	25 / case
364.057	40	5	1.7	1.0	0.39	Grey	25 / case
364.067	40	6	2.0	1.2	0.54	Light green	25 / case
364.087	40	8	2.5	1.5	0.81	Blue	25 / case
361.047	50	4	1.5	0.8	0.32	Purple	50 / case
361.057	50	5	1.7	1.0	0.47	Grey	50 / case
361.067	50	6	2.0	1.2	0.65	Light green	50 / case
361.087	50	8	2.5	1.5	0.99	Blue	50 / case
361.107	50	10	3.3	2.0	1.7	Black	50 / case
363.047	75	4	1.5	0.8	0.45	Purple	25 / case
363.057	75	5	1.7	1.0	0.67	Grey	25 / case
363.067	75	6	2.0	1.2	0.94	Light green	25 / case
363.087	75	8	2.5	1.5	1.4	Blue	25 / case
363.107	75	10	3.3	2.0	2.5	Black	25 / case
363.127	75	12	4.0	2.5	3.9	White	25 / case
366.087	90	8	2.5	1.5	1.7	Blue	25 / case
366.107	90	10	3.3	2.0	3.0	Black	25 / case
366.127	90	12	4.0	2.5	4.6	White	25 / case
366.147	90	14	4.5	3.0	6.6	Dark green	25 / case
366.167	90	16	5.0	3.5	8.9	Orange	25 / case
366.187	90	18	6.0	4.0	11.7	Red	25 / case
362.047	125	4	1.5	0.8	0.7	Purple	25 / case
362.057	125	5	1.7	1.0	1.06	Grey	25 / case
362.067	125	6	2.0	1.2	1.5	Light green	25 / case
362.087	125	8	2.5	1.5	2.3	Blue	25 / case
362.107	125	10	3.3	2.0	4.0	Black	25 / case
362.127	125	12	4.0	2.5	6.3	White	25 / case
362.147	125	14	4.5	3.0	9.0	Dark green	25 / case
362.167	125	16	5.0	3.5	12.3	Orange	25 / case
362.187	125	18	6.0	4.0	16.0	Red	25 / case

## POLYURETHANE FEEDING TUBES WITHOUT STYLET

Adapted for children and adults, the PUR nutritif feeding tubes are intended to be used for medium-term (1 to 4 weeks) and feature the following characteristics:

- Made of PUR
- Closed distal tip with 2 lateral eyes (4Fr to 12Fr) and 4 lateral eyes (14Fr to 18Fr)
- Small cap with color code to quickly identify the size
- Numerical markings every centimeter:
  - from 5cm to 35cm for codes 1364, 1361 & 1363
  - from 5cm to 70cm for codes 1366, 1362 & 1365
- Radiopaque line from 4Fr to 6Fr; and 3/4 radiopaque tubing from 8Fr
- Transparent to better visualize the administered solution.




Code	PUR tube without stylet					Color code	Quantity
	Length cm	Fr	Ext. Ø mm	Int. Ø mm	Prim. vol. ml		Units/case
1364.047	40	4	1.5	0.8	0.27	Purple	25 / case
1364.057	40	5	1.7	1.0	0.39	Grey	25 / case
1364.067	40	6	2.0	1.2	0.54	Light green	25 / case
1364.087	40	8	2.5	1.5	0.81	Blue	25 / case
1361.047	50	4	1.5	0.8	0.32	Purple	50 / case
1361.057	50	5	1.7	1.0	0.47	Grey	50 / case
1361.067	50	6	2.0	1.2	0.65	Light green	50 / case
1361.087	50	8	2.5	1.5	0.99	Blue	50 / case
1361.107	50	10	3.3	2.0	1.7	Black	50 / case
1363.047	75	4	1.5	0.8	0.45	Purple	25 / case
1363.057	75	5	1.7	1.0	0.67	Grey	25 / case
1363.067	75	6	2.0	1.2	0.94	Light green	25 / case
1363.087	75	8	2.5	1.5	1.4	Blue	25 / case
1363.107	75	10	3.3	2.0	2.5	Black	25 / case
1363.127	75	12	4.0	2.5	3.9	White	25 / case
1366.087	90	8	2.5	1.5	1.7	Blue	25 / case
1366.107	90	10	3.3	2.0	3.0	Black	25 / case
1366.127	90	12	4.0	2.5	4.6	White	25 / case
1362.047	125	4	1.5	0.8	0.7	Purple	25 / case
1362.057	125	5	1.7	1.0	1.06	Grey	25 / case
1362.067	125	6	2.0	1.2	1.5	Light green	25 / case
1362.087	125	8	2.5	1.5	2.3	Blue	25 / case
1362.107	125	10	3.3	2.0	4.0	Black	25 / case
1362.127	125	12	4.0	2.5	6.3	White	25 / case
1362.147	125	14	4.5	3.0	9.0	Dark green	45 / case
1362.167	125	16	5.0	3.5	12.3	Orange	45 / case

## POLYURETHANE FEEDING TUBES WITH STYLET

Adapted for children and adults, the PUR nutritif feeding tubes are intended to be used for medium-term enteral feeding (1 to 4 weeks) and feature the following characteristics:

- Made of PUR
- Stylet protected by a ball end to facilitate the placement
- Distal tip with 2 lateral eyes (6Fr to 12Fr) and 4 lateral eyes (14Fr to 18Fr)
- Small cap with color code to quickly identify the size
- Numerical markings every centimeter
  - from 5cm to 35cm for code 1363
  - from 5cm to 70cm for codes 1366 & 1362
- Radiopaque.



Code	PUR tube with stylet						Color code	Quantity
	Length cm	Fr	Ext. Ø mm	Int. Ø mm	Distal tip	Prim. vol. ml		Units/case
1366.0671	90	6	2.0	1.2	Closed	1.1	Light green	25 / case
1366.0871	90	8	2.5	1.7	Closed	2.2	Blue	25 / case
1366.1071	90	10	3.3	2.5	Closed	4.6	Black	25 / case
1366.1271	90	12	4.0	3.2	Closed	7.5	White	25 / case
1366.1471	90	14	4.5	3.5	Closed	9.0	Dark green	25 / case
1366.1671	90	16	5.3	4.0	Closed	11.7	Orange	20 / case
1366.1871	90	18	6.0	4.5	Closed	14.8	Red	20 / case
1366.0872	90	8	2.5	1.7	Open	2.2	Blue	25 / case
1366.1072	90	10	3.3	2.5	Open	4.6	Black	25 / case
1366.1272	90	12	4.0	3.2	Open	7.5	White	25 / case
1362.0671	125	6	2.0	1.2	Closed	1.5	Light green	25 / case
1362.0871	125	8	2.5	1.7	Closed	2.9	Blue	25 / case
1362.1071	125	10	3.3	2.5	Closed	6.3	Black	25 / case
1362.1271	125	12	4.0	3.2	Closed	10.3	White	25 / case
1362.1471	125	14	4.5	3.5	Closed	12.3	Dark green	25 / case
1362.1671	125	16	5.3	4.0	Closed	16.0	Orange	20 / case
1362.1871	125	18	6.0	4.5	Closed	20.3	Red	20 / case
1362.0872	125	8	2.5	1.7	Open	2.9	Blue	25 / case
1362.1072	125	10	3.3	2.5	Open	6.3	Black	25 / case

## SYRINGES

### LOW DOSE TIP SYRINGES

The Low Dose Tip (LDT) **nutrifit** syringes incorporate a purple-coloured plunger and "Enteral" print, enabling instant identification for enteral use. They are sterile and individually packed.



Code	Transparent LDT syringes			Cap	Quantity
	Volume	Graduation	ENFit® connection		Units/case
1015.018	1 ml	0.01 ml	Centered	No	100 / case
1015.028	2.5 ml	0.1 ml	Centered	No	100 / case
1015.058	5 ml	0.2 ml	Centered	No	100 / case

## SYRINGES

### TRANSPARENT SYRINGES

The transparent **nutrifit** syringes incorporate a purple-coloured plunger and "Enteral" print, enabling instant identification for enteral use. They are sterile and individually packed.



Code	Transparent syringes			Cap	Quantity
	Volume	Graduation	ENFit® connection		Units/case
1015.107	10 ml	0.5 ml	Offset	No	100 / case
1015.207	20 ml	0.5 ml	Offset	No	80 / case
1015.607	60 ml	1 ml	Offset	No	60 / case
1015.1007	100 ml	2 ml	Centered	No	50 / case

## SYRINGES ACCESSORIES

### MALE CAPS FOR **nutrifit** SYRINGES

These caps are specifically designed for **nutrifit** syringes. Purple-coloured, it enables instant identification for enteral use. Compatible with Low Dose Tip syringes.



Code	Packaging	Units/case
828.017	Individually packed	200 / case
828.087	Blister containing 8 caps	200 / case

### VENTED CAPS

Vented cap are used for **nutrifit** syringes (20 and 60ml) to enable gravity feeding and protecting the solution of enteral feed from contamination.

Code 818.07  
Quantity Case: 30




## ACCESSORIES

### EXTENSION TUBES

#### STANDARD EXTENSION TUBES

These nutri**fit** extension tubes keep the feeding tube away from the syringe pump or from the feeding pump to improve the patient's comfort.

- Purple-coloured PVC extension tube
- Attached cap to close the tube and avoid contamination depending on the reference.




Code	Tube					Quantity
	Length cm	Ext. Ø mm	Int. Ø mm	Priming volume ml	Cap	Units/case
368.097	90	2.5	1.5	1.38	No	200 / case
368.107	100	2.5	1.5	1.9	No	45 / case
368.157	150	2.5	1.5	2.8	No	45 / case
368.1572	150	2.5	1.5	2.3	No	200 / case

#### EXTENSION TUBES WITH CLAMP

This nutri**fit** device gives extended connections to pumps and allows greater freedom of movement for ambulatory patients:

- Purple - coloured PVC extension tube
- With Roberts clamp
- Attached cap to close the tube and avoid contamination.

Code	Tube				Quantity
	Length cm	Ext. Ø mm	Int. Ø mm	Priming volume ml	Units/case
368.1071	100	2.5	1.5	1.9	45 / case
368.1571	150	2.5	1.5	2.8	45 / case



## SAMPLING DEVICES

### DRAWING-UP STRAWS

The nutri**fit** drawing-up straw is a device connected to a nutri**fit** syringe for medication sampling.



Code	Drawing-up straw				Quantity
	Length cm	Ext. Ø mm	Int. Ø mm	Priming volume ml	Units/case
817.057	5	2.0	1.2	0.13	100 / case
817.0572	5	2.0	1.0	0.10	200 / case
817.157	15	3.3	2.0	0.57	50 / case
817.207	20	4.1	2.8	1.3	200 / case
817.208	20	7.4	4.9	3.75	200 / case

### BOTTLE CAPS FOR SAMPLING

These nutri**fit** caps secure the medicine bottle and allow the sampling with a nutri**fit** syringe.



Code	Bottle caps for sampling	Quantity
	Ø mm	Units/case
820.1297	12.9	30 / case
820.2047	20.4	30 / case
820.2777	27.7	30 / case

## OTHER ACCESSORIES

### CHRISTMAS TREE ADAPTERS

These nutri**fit** adapters help customize your sets and help manage the ENFit® transition period.



Code	Product	Quantity Units/case
0VEACC406	Male ENFit® / Christmas tree adaptor (without cap, non-sterile)	50 / case
0VEACC416	Male ENFit® / Christmas tree adaptor (with cap, non-sterile)	50 / case