


# Esophageal Balloon Catheter Set

## For Intrathoracic Pressure Monitoring

- Ease of Insertion - The catheter comes with a new, easy-to-remove PTFE-coated stylet to aid the insertion procedure
- Comfort - The 5 French catheter is soft and pliable, offering maximum patient comfort over extended time periods
- Ease of Positioning - The new proximal Y-connector allows for precise balloon positioning with the stylet in place
- Depth Markings - The catheter is radio opaque with numerical depth markings to facilitate positioning
- Ease of Connection - The catheter has a standard female Luer connector
- Several Applications - Ideal for monitoring ventilator and sleep apnea patients and calculating lung compliance



Item #	Description
47-9005	5 Fr Catheter, latex-free Balloon, Stylet, 36" Pressure Monitoring Extension, 3-way Stopcock (5/Box)

