

# exchange transfusion tray

Set containing the accessories necessary to perform removal and replacement of a patient's blood.



Product Code	Contents	Units/Case
275.001	1x four-way stopcock equipped with 1 injection site (for injections of additional medication)	1
	1x extension tube for evacuation of the discarded blood	
	1x 5Fr exchange transfusion catheter (PVC - XRO) 37 cm long, with 2 sides holes and centimetric markings	
	1x 7Fr exchange transfusion catheter (PVC - XRO) 37 cm long, with 2 sides holes and centimetric markings	
	2x 20 ml Luer-lock syringes	
	1x 10 ml Luer-lock syringe	
	1x 16x0.5 mm (25G) Luer hypodermic needle	
	1x graduated plastic container (2000 ml)	
	1x transfusion set, 163 cm long, 18.8 ml dead volume	
	1x 15 cm ruler for venous pressure measurement	
	3x 50 x 50 mm gauze dressings	
	1x 50 x 60 cm fenestrated drape	
	1x air-venting needle	
	1x control sheet	