

# gastro-duodenal

## Feeding Tubes

125 cm tubes for placement in the duodenum.

Weighted tubes: Silicone (for both medium and long-term use) tube with weighted metal tip to facilitate placement

Levin type tubes: DEHP-free PVC (for short-term use) tubes for gastric decompression or gavage feeding.



**Silicone tube**



**Levin type tubes with standard connector**

### Weighted tubes (silicone)

Product Code	Size	Units/Case
2395.05	5 Fr	25
2395.06	6 Fr	25
2395.09	9 Fr	25

### Levin type tubes with standard connector (PVC)

Product Code	Size	Units/Case
391.10	10 Fr	50
391.12	12 Fr	50
391.14	14 Fr	50
391.16	16 Fr	50