

H/S Catheters


For Enhanced Visualization of Uterine Pathology

H/S CATHETER SET

These easy-to-use devices provide enhanced visualization. Designed for use in sonohysterography and hysterosalpingography, the catheter can be used with oil- or saline-based media depending on the procedure requirements.

Features and Benefits

- Available with integrated stylet for easy insertion
- Unique placement sheath eliminates the need for a tenaculum
- Comes in both 5 and 7 French
- Ideal for office use
- Latex-Free construction



Pre-Saline Infusion

Post-Saline Infusion

Item #	Description
61-5005	H/S Catheter, Kraton, 5 Fr (10/Box)
61-5007	H/S Catheter, Kraton, 7 Fr (10/Box)
61-0020	H/S Stylet Catheter Stiffener (10/Box)
61-3005	H/S Catheter with integrated Stylet 5FR (10/Box)
61-3007	H/S Catheter with integrated Stylet 7FR (10/Box)