

# INCA<sup>®</sup>

## Infant Nasal CPAP Assembly

INCA<sup>®</sup> Infant Nasal CPAP Assembly is the most respected Continuous Positive Airway Pressure (CPAP) device among clinicians. INCA affords all the benefits of CPAP/Bubble CPAP for newborns in your care.

Features enhanced precision fit for neonates:

- Five anatomically proportional nasal cannula sizes for neonates
- Clear, extra-soft tubing provides custom shaping and maximal stability
- Chin strap is available to help prevent air leakage around the mouth
- Available in Complete and Replacement Sets



### INCA COMPLETE SET

Consists of a nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters, one in-line "T" pressure adapter, a knit cap with hook/loop strap, two foam blocks and a hook/loop chin strap. One INCA Sizing Gauge is included in every box of five sets.



Item #	Description
44-2707	Set with 7.5 Fr Cannula (5/Box)
44-2709	Set with 9 Fr Cannula (5/Box)
44-2710	Set with 10.5 Fr Cannula (5/Box)
44-2712	Set with 12 Fr Cannula (5/Box)
44-2715	Set with 15 Fr Cannula (5/Box)

### INCA REPLACEMENT SET

Includes the nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters and one in-line pressure "T" adapter.




Item #	Description
44-0707	Set with 7.5 Fr Cannula (5/Box)
44-0709	Set with 9 Fr Cannula (5/Box)
44-0710	Set with 10.5 Fr Cannula (5/Box)
44-0712	Set with 12 Fr Cannula (5/Box)
44-0715	Set with 15 Fr Cannula (5/Box)

### INCA ACCESSORY PACK

Includes five knit caps with hook/loop strap, five chin straps with ten hook fasteners (not attached) and ten foam blocks.

1 Pack = 5 Hats and 10 Foam Blocks



Item #	Description
44-2001	X-Small Pack (5 XS Caps, 10 Foam Blocks)
44-2002	Small Pack (5 S Caps, 10 Foam Blocks)
44-2003	Medium Pack (5 M Caps, 10 Foam Blocks)
44-2004	Large Pack (5 L Caps, 10 Foam Blocks)

### INCA NASAL SIZING GAUGE

For appropriate nasal prong fitting, use the INCA Sizing Gauge. Wipe gauge with sterile 70 percent isopropyl alcohol prior to each use. Hold the gauge up to the nasal area to determine appropriate cannula size.



Item #	Description
025-001	INCA Nasal Sizing Gauge (1/Pack)