

Rapid Infusion

INFU-SURG® CLARITY PRESSURE INFUSION BAG

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ethox



Clear panel enables full visibility of solution bag to monitor fluid level and easily scan patient barcode



- Transparent panel for clear visibility of solution bag
- Ideal for intravenous and intra-arterial infusion treatment including A-line pressure monitoring
- Dedicated hook enables quick change of solution bag while providing proper fluid bag positioning in pocket
- Pressure relief valve prevents over-inflation (330 mm Hg)
- Color-coded gauge provides accurate pressure monitoring (150-300 mm Hg)
- Single-handed inflation and deflation with oval-shaped bulb makes it easy to use and requires minimal training
- Welded panel and hook construction are more durable than sewn designs to withstand repeated use
- Reinforced hanging loop allows for quick and easy hanging on IV pole
- Three-way stopcock ensures precise control of pressure
- Piston-style gauge allows for 360° visibility



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation

ITEM	DESCRIPTION	VOLUME	PK
4050	Infu-Surg® Clarity	500 mL	5
4100	Infu-Surg® Clarity	1000 mL	25