



Rapid Infusion

INFU-SURG® PRESSURE INFUSION BAG

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ethox



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation

- Ideal for intravenous and intra-arterial infusion treatment including A-line pressure monitoring
- Pressure relief valve prevents over inflation (330 mm Hg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- Dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole
- Reinforced hanging loop allows for quick and easy hanging on IV pole
- Color-coded gauge makes for accurate pressure monitoring (150-300 mm Hg)
- Three-way stopcock ensures precise control of pressure
- Piston-style gauge allows for 360° visibility

ITEM#	DESCRIPTION	VOLUME	PK
4005H	Infu-Surg® with Hook	500 mL	5
4010H	Infu-Surg® with Hook	1000 mL	5
4030	Infu-Surg®	3000 mL	5



DISTRIBUTED BY:
ADVANCED SURGI-PHARM INC.
www.advancedsurgipharma.com
Tel:1-800-661-5432