

IV Accessories

Ramp of stopcocks (1 way)

Product Code	Description	Extension tube			Lateral accesses	Units/Box
		Ext. Ø	Int. Ø	Length	number	
873.24	Ramp of stopcocks (1 way)	4.0 mm	2.5 mm	200 cm	4	15



Ramp of Stopcocks (1 way): female Luer-lock/female Luer-lock

Stopcock

Product Code	Description	Units/Box
872.10	Lockable one-way stopcock	50



Lockable one-way stopcock

Connectors

Product Code	Description	Units/Box
892.00	Female Luer-lock to female Luer-lock connector	30
893.00	Male Luer Lock to Male Luer Lock Connector	30



Female Luer Lock Connector



Male Luer-lock to male Luer-lock connector

Extension tubes (single lumen)

Product Code	Description	Extension tube			Units/Box
		Ext. Ø	Int. Ø	Length	
220.050	Lectroflex extension tube	4.0 mm	2.5 mm	50 cm	25
221.100	Male Luer-lock to male Luer-lock extension tube	4.0 mm	2.5 mm	100 cm	25



Lectroflex



Male Luer-lock to male Luer-lock extension tube