

## K-Wires & Steinmann Pins

Packaging: 6 or 10 units per plastic tube / 1 tube per package

### K-Wires

	DOUBLE TROCAR	SINGLE TROCAR	Length (in)	Length (cm)	Diameter (in)	Diameter (mm)
SMOOTH	Item #	Item #				
	SU78-2000	-	4"	10 cm	.031"	0.8 mm
	SU78-2010	-	4"	10 cm	.039"	1.0 mm
	SU78-2020	-	4"	10 cm	.047"	1.2 mm
	SU78-2040	-	4"	10 cm	.054"	1.4 mm
	SU78-2030	-	4"	10 cm	.062"	1.6 mm
	SU78-1210	-	6"	15 cm	.031"	0.8 mm
	SU78-1220	SU78-2810	6"	15cm	.039"	1.0 mm
	SU78-1230	SU78-2820	6"	15 cm	.047"	1.2 mm
	SU78-1240	-	6"	15 cm	.054"	1.4 mm
	SU78-1250	SU78-2850	6"	15 cm	.062"	1.6 mm
	SU78-2105	SU78-2500	9"	23 cm	.031"	0.8 mm
	SU78-2110	SU78-2510	9"	23 cm	.039"	1.0 mm
	SU78-2120	SU78-2520	9"	23 cm	.047"	1.2 mm
	SU78-2140	-	9"	23 cm	.054"	1.4 mm
SU78-2130	SU78-2530	9"	23 cm	.062"	1.6 mm	
FULL THREAD	SU78-4030	-	9"	23 cm	.062"	1.6 mm

### DOUBLE DIAMOND SINGLE DIAMOND

	DOUBLE DIAMOND	SINGLE DIAMOND	Length (in)	Length (cm)	Diameter (in)	Diameter (mm)
SMOOTH	Item #	Item #				
	-	SU78-2810	6"	15 cm	.039"	1.0 mm
	-	SU78-2820	6"	15 cm	.047"	1.2 mm
	-	SU78-2850	6"	15 cm	.062"	1.6 mm
	SU78-3110	SU78-3410	9"	23 cm	.039"	1.0 mm
	SU78-3120	SU78-3420	9"	23 cm	.047"	1.2 mm
SU78-3130	SU78-3430	9"	23 cm	.062"	1.6 mm	
FULL THREAD	-	SU78-4070	9"	23 cm	.062"	1.6 mm

## K-Wires & Steinmann Pins

### Steinmann Pins

**DOUBLE TROCAR      SINGLE TROCAR**

	Item #	Item #	Length (in)	Length (cm)	Diameter (in)	Diameter (mm)
<b>SMOOTH</b>	-	SU78.6525	6"	15 cm	.098"	2.5 mm
	-	SU78.6429	6"	15 cm	.110"	2.8 mm
	SU78.5500	SU78.6100	9"	23 cm	.078"	2.0 mm
	SU78.5530	SU78.6130	9"	23 cm	.098"	2.5 mm
	SU78.5560	SU78.6160	9"	23 cm	.110"	2.8 mm
	SU78.5590	SU78.6190	9"	23 cm	.118"	3.0 mm
	SU78.5620	SU78.6220	9"	23 cm	.141"	3.5 mm
	SU78.5650	SU78.6250	9"	23 cm	.157"	4.0 mm
	SU78.5680	SU78.6280	9"	23 cm	.177"	4.5 mm
	SU78.5690	SU78.6290	9"	23 cm	.197"	5.0 mm
<b>FULL THREAD</b>	SU78.8500	-	9"	23 cm	.078"	2.0 mm
	SU78.8530	-	9"	23 cm	.098"	2.5 mm
	SU78.8560	-	9"	23 cm	.110"	2.8 mm
	SU78.8590	-	9"	23 cm	.118"	3.0 mm
	SU78.8620	-	9"	23 cm	.141"	3.5 mm
	SU78.8650	-	9"	23 cm	.157"	4.0 mm

**DOUBLE DIAMOND      SINGLE DIAMOND**

	Item #	Item #	Length (in)	Length (cm)	Diameter (in)	Diameter (mm)
<b>SMOOTH</b>	SU78.7000	SU78.7780	9"	23 cm	.078"	2.0 mm
	SU78.7060	SU78.7660	9"	23 cm	.110"	2.8 mm
	-	SU78.7690	9"	23 cm	.118"	3.0 mm
	-	SU78.7720	9"	23 cm	.141"	3.5 mm
	-	SU78.7782	9"	23 cm	.177"	4.5 mm
	SU78.1350	-	6"	15 cm	.141"	3.5 mm
<b>FULL THREAD</b>	SU78.8300	-	9"	23 cm	.078"	2.0 mm
	-	SU78.8030	9"	23 cm	.098"	2.5 mm
	-	SU78.8120	9"	23 cm	.141"	3.5 mm
<b>3 SHANK END</b>	-	SU78.7630	9"	23 cm	.110"	2.8 mm

## K-Wires & Steinmann Pins

	DOUBLE TROCAR	SINGLE TROCAR
Smooth		
Partial Thread		
Full Thread		
Smooth - Round Hole		
	DOUBLE DIAMOND	SINGLE DIAMOND
Smooth		
Partial Thread		
Full Thread		
Smooth - Round Hole		
3 Shank End		

## Sterilization Racks

Sterilization racks are designed for use during washing, sterilizing, storing/transporting, and dispensing K-Wires, Steinmann Pins, or Twist Drills. The racks are lightweight but made to withstand repeated use. They fold closed for sterilization and easily open to convert to a table top stand for use in the O.R. Dimensions are clearly marked on the rack to assist with quick identification.



### K-Wire and Pin Rack

Number	Details	Dimensions	Quantity/ Package
SU98.5409-RO	Anodized aluminum, holds 9" K-Wires and pins Holds 6 of each of the following diameters: 0.028" (0.7 mm) 0.045" (1.1 mm) 0.062" (1.6 mm) 0.093" (2.4 mm) 0.125" (3.2 mm) 0.156" (4.0 mm) 0.035" (0.9 mm) 0.054" (1.4 mm) 0.078" (2.0 mm) 0.109" (2.8 mm) 0.140" (3.5 mm) 0.177" (4.5 mm)	12.5" x 5.25" x 1.5" (32 x 13 x 3.8)	1



### K-Wire Rack

Number	Details	Dimensions	Quantity/ Package
SU98.5404-RO	Anodized aluminum, holds 4-6" K-Wires and pins Holds 6 of each of the following diameters: 0.028" (0.7 mm) 0.045" (1.1 mm) 0.062" (1.6 mm) 0.035" (0.9 mm) 0.054" (1.4 mm)	12.5" x 5.25" x 1.5" (32 x 13 x 3.8)	1