

# leadercath


## Arterial LeaderCath (PE)

Although it can be inserted in a variety of arteries, the most common sites for the placement of this catheter are the radial and the femoral artery.

The Seldinger technique makes arterial puncture easier particularly in hypotensive patients with a weak peripheral pulse and when arteries are atherosclerosed or narrow.

The product is supplied in a rigid blister pack containing:

- 1 introducer needle
- 1 straight guidewire
- 1 transparent and radiopaque catheter in polyethylene with winged hub



Product Code	Catheter						Puncture Needle				Guidewire		Unit/Box
	Int. Ø mm	Ext. Ø mm	Length cm	Gauge G	Ø Fr	Flow rate ml/min	Int. Ø mm	Ext. Ø	Gauge G	Length mm	Ø mm	Length mm	
115.090	0.6	0.9	8	20	3	24	0.6	0.9	20	38	0.53	20	20
115.11	0.8	1.2	10	18	4	24	0.8	1.1	19	54	0.71	30	20