

# lectrocath

A range of polyethylene extension lines, suitable for arterial pressure monitoring and venous medication delivery.

The advantages include:

- a Polyethylene material which allows for compatibility with commonly used drugs and reduces adsorption, absorption or permeation that may occur with PVC
- Latex-Free
- DEHP-Free
- Pyrogen-Free



Feature	Benefit
PE tubing	Drug-compatible
Colour-coded clips	For arterial or venous identification

Product Code	ID - OD	Length	Capacity	Units/Box
1155.01	1.0 - 2.0 mm	11 cm	0.4 ml	50
1155.03	1.0 - 2.0 mm	30 cm	0.6 ml	50
1155.05	1.0 - 2.0 mm	50 cm	0.7 ml	50
1155.10	1.0 - 2.0 mm	100 cm	1.2 ml	50
1155.15	1.0 - 2.0 mm	150 cm	1.6 ml	50
1155.20	1.0 - 2.0 mm	200 cm	2.0 ml	50
1159.01	1.5 - 2.5 mm	11 cm	0.5 ml	50
1159.03	1.5 - 2.5 mm	30 cm	0.8 ml	50
1159.05	1.5 - 2.5 mm	50 cm	1.1 ml	50
1159.10	1.5 - 2.5 mm	100 cm	2.0 ml	50
1159.15	1.5 - 2.5 mm	150 cm	2.8 ml	50
1159.20	1.5 - 2.5 mm	200 cm	3.6 ml	50