

lectrospiral[®]

Coiled extension tube

Lectrospiral[®] is a flexible coiled extension tube that provides patients with freedom of movement.

Flexible

- Coiled tubing expands and contracts providing patients with greater freedom of movement than straight extension lines.
- Ideal to use with pediatric, home care and oncology patients.

Safe

- Polyethylene material allows for compatibility with commonly infused drugs and reduces adsorption, absorption or permeation that may occur with PVC.
- Extension tubing minimizes potential of trauma at IV site.



Product Code	Length	Inner Diameter	Outer Diameter	Priming Volume	Units/Case
I 155.60	100 cm	1.0 mm	2.0 mm	1.1 ml	15
I 155.65	150 cm	1.0 mm	2.0 mm	1.5 ml	15
I 155.70	200 cm	1.0 mm	2.0 mm	2.0 ml	15
I 155.80	300 cm	1.0 mm	2.0 mm	2.8 ml	15
I 155.90	400 cm	1.0 mm	2.0 mm	3.7 ml	15
I 159.60	100 cm	1.5 mm	2.5 mm	1.9 ml	15
I 159.65	150 cm	1.5 mm	2.5 mm	2.7 ml	15
I 159.70	200 cm	1.5 mm	2.5 mm	3.6 ml	15
I 159.80	300 cm	1.5 mm	2.5 mm	5.2 ml	15
I 159.90	400 cm	1.5 mm	2.5 mm	6.9 ml	15