

STERILE SURGICAL LAMINATED ULTRASONICALLY WELDED GOWNS

LEVEL 3
STERILE
ULTRASONICALLY WELDED

Our **Level 3 Disposable Sterile Surgical Laminated Gowns** provide high performance with a fluid resistant barrier protection.

- Fabric Compliant to **EN 20811** and **EN ISO 22610**
- Tested and compliant to **AAMI PB70 - Level 3 Liquid Barrier Performance**
- Chemo Resistance to **12 CHEMOTHERAPY DRUGS: BICNU, Adriamycin, Cytosan, Taxol, Toposar, 5-Flu, Thiotepa, Cisplatin, DTIC, Mitomycin C, Oncovin, and Amethopterin**
- 56 gr/m² fabric for maximum **protection, breathability, and comfort**
- Knitted soft wrist cuffs
- Tie closure
- **Ultrasonically welded seams**
- **Sterile and aseptically packed**

Individually packed, each gown comes in a set designed to be used during surgery to **reduce contamination and infection** between the patient and medical team.

Each set includes:

- 1 x Sterile Surgical Laminated Gown - Level 3 (Size S - XXXL)
- 1 x Blue Sterilization Wrap (60 cm x 60 cm / 16" x 16")
- 1 x White Cellulose Hand Towel (30 cm x 40 cm / 12" x 16")

PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
41.001.01	S	110 cm / 43"	55 cm / 22"	140 cm / 55"
41.002.01	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
41.003.01	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
41.004.01	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
41.005.01	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"
41.006.01	XXXL	145 cm / 57"	66 cm / 26"	160 cm / 63"



PACKAGING: 40 INDIVIDUALLY BAGGED SETS / BOX