

BATRIK

LEVEL 4 PROTECTIVE LAMINATED ULTRASONICALLY WELDED & REINFORCED GOWNS

LEVEL 4
NON-STERILE
ULTRASONICALLY WELDED
& REINFORCED SEAMS

Our **Level 4 Disposable Protective Laminated Gowns** provide high performance for use in high risk situations and/or when infectious diseases are suspected.

- Fabric Compliant to **EN 20811** and **EN ISO 22610**
- Tested and compliant to **AAMI PB70 - Level 4 Liquid Barrier Performance and blood and viral penetration resistance**
- 56 gr/m² soft non-woven spunbond fabric coated with a polyethylene film provides maximum **protection, breathability, and comfort**
- Tested and compliant to **ASTM F739-12**
- Knitted soft wrist cuffs
- Tie closure
- **Ultrasonically welded & reinforced seams**
- **Non-Sterile**



**INSIDE VIEW OF
ULTRASONICALLY WELDED &
REINFORCED SEAMS**

PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
SF-402-S	S	110 cm / 43"	55 cm / 22"	140 cm / 55"
SF-402-M	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
SF-402-L	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
SF-402-XL	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
SF-402-XXL	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"
SF-402-XXXL	XXXL	145 cm / 57"	66 cm / 26"	160 cm / 63"

PACKAGING: 50 GOWNS / BOX