

# LEVEL 4 PROTECTIVE LAMINATED ULTRASONICALLY WELDED & REINFORCED GOWNS

**LEVEL 4**  
**NON-STERILE**  
 ULTRASONICALLY WELDED  
 & REINFORCED SEAMS

Our **Level 4 Disposable Protective Laminated Gowns** provide high performance for use in high risk situations and/or when infectious diseases are suspected.

- Fabric Compliant to **EN 20811** and **EN ISO 22610**
- Tested and compliant to **AAMI PB70 - Level 4 Liquid Barrier Performance** and **blood and viral penetration resistance**
- 56 gr/m<sup>2</sup> soft non-woven spunbond fabric coated with a polyethylene film provides maximum **protection, breathability, and comfort**
- Tested and compliant to **ASTM F739-12**
- Knitted soft wrist cuffs
- Tie closure
- **Ultrasonically welded & reinforced seams**
- **Non-Sterile**



INSIDE VIEW OF  
ULTRASONICALLY WELDED &  
REINFORCED SEAMS

## PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
<b>SF-402-S</b>	<b>S</b>	110 cm / 43"	55 cm / 22"	140 cm / 55"
<b>SF-402-M</b>	<b>M</b>	115 cm / 45"	58 cm / 23"	145 cm / 57"
<b>SF-402-L</b>	<b>L</b>	127 cm / 50"	60 cm / 24"	150 cm / 59"
<b>SF-402-XL</b>	<b>XL</b>	135 cm / 53"	62 cm / 24"	150 cm / 59"
<b>SF-402-XXL</b>	<b>XXL</b>	140 cm / 55"	64 cm / 25"	155 cm / 61"
<b>SF-402-XXXL</b>	<b>XXXL</b>	145 cm / 57"	66 cm / 26"	160 cm / 63"

PACKAGING: 50 GOWNS / BOX