

# STERILE SURGICAL LAMINATED ULTRASONICALLY WELDED & REINFORCED GOWNS

**LEVEL 4  
STERILE**  
ULTRASONICALLY WELDED  
& REINFORCED SEAMS

Our **Level 4 Disposable Sterile Surgical Laminated Gowns** provide high performance with a fluid resistant barrier protection.

- Fabric Compliant to **EN 20811** and **EN ISO 22610**
- Tested and compliant to **AAMI PB70 - Level 4 Liquid Barrier Performance** and **blood and viral penetration resistance**
- 56 gr/m<sup>2</sup> soft non-woven spunbond fabric coated with a polyethylene film provides maximum **protection, breathability, and comfort**
- Tested and compliant to **ASTM F739-12**
- Knitted soft wrist cuffs
- Tie closure
- **Ultrasonically welded & reinforced seams**
- **Sterile and aseptically packed**

Individually packed, each gown comes in a set designed to be used during surgery to **reduce contamination and infection** between the patient and medical team.

## Each set includes:

- 1 x Sterile Surgical Laminated Gown - Level 4 (Size S - XXXL)
- 1 x Blue Sterilization Wrap (60 cm x 60 cm / 16" x 16")
- 1 x White Cellulose Hand Towel (30 cm x 40 cm / 12" x 16")



**INSIDE VIEW OF  
ULTRASONICALLY  
WELDED &  
REINFORCED SEAMS**

## PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
<b>46.001.01</b>	<b>S</b>	110 cm / 43"	55 cm / 22"	140 cm / 55"
<b>46.002.01</b>	<b>M</b>	115 cm / 45"	58 cm / 23"	145 cm / 57"
<b>46.003.01</b>	<b>L</b>	127 cm / 50"	60 cm / 24"	150 cm / 59"
<b>46.004.01</b>	<b>XL</b>	135 cm / 53"	62 cm / 24"	150 cm / 59"
<b>46.005.01</b>	<b>XXL</b>	140 cm / 55"	64 cm / 25"	155 cm / 61"
<b>46.006.01</b>	<b>XXXL</b>	145 cm / 57"	66 cm / 26"	160 cm / 63"

**PACKAGING: 40 INDIVIDUALLY BAGGED SETS / BOX**