

BATRIK

LILA'S INSTRUMENT STAIN REMOVER

Lila's special blend of naturally sourced powders provides outstanding removal of stubborn stains.

Convenient-to-use formula, you can achieve spotless results every time!

It can be used on surgical/dental instruments, autoclaves, ultrasonic cleaners, stainless steel countertops, trays, bowls and more.

Specifically designed for use on stainless steel, chrome, nickel, brass, copper, and silver.

**LSR-8501**


EASE TO USE AS 1,2,3...

1. Dampen a cloth or grab a Lila's Low Linting Dura-Towel and dab into the jar.
2. Gently rub the area that is stained.¹
3. Rinse thoroughly for a flawless shine.²

A little Goes a long Way!

Even great to tackle stubborn stains on steel pots and pans in your kitchen with ease.

Lila's - If you have used Spectra-Scrub or Integra's instrument stain remover then you will love Lila's!

**60-46003****17-DT1213C**

LILA'S SCRUB PAD

- Say goodbye to extra stubborn stains with Lila's Scrub Pad without leaving any metal residue.
- Simply add a small amount of Lila's Stain Remover with one of our Scrub Pad's for very large areas.
- The gentle yet effective scrub is the perfect choice for Ophthalmic Instruments.¹

LILA'S LOW-LINTING DURA-TOWEL

- Lila's Dura-Towel has long lasting strength and has a minimal lint composition.
- Its ideal dimensions and unique attributes guarantee optimal performance, setting it apart from ordinary paper towels.

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
LSR-8501	LILA'S INSTRUMENT STAIN REMOVER 3 oz (85g) plastic jar with a convenient screw-on lid	Ø 2.8" / 70 mm	1 Jar / Bag
LSR-8512	LILA'S INSTRUMENT STAIN REMOVER 3 oz (85g) plastic jar with a convenient screw-on lid	Ø 2.8" / 70 mm	12 Jars / Box
17-DT1213C	LILA'S LOW LINTING DURA-TOWEL 1080 towels	12" x 13" / 31 cm x 33 cm Folded: 3" x 3.25" / 7.5 cm x 8.25 cm	90 / Bag 12 Bags / Box
60-46003	LILA'S SCRUB PAD	4" x 6" / 10 cm x 15 cm	3 / Bag
LSR-BUNDLE	LILA'S STAIN REMOVER + LILA'S DURA-TOWEL 3 oz (85g) plastic jar + 90 towels	Ø 2.8" / 70 mm 12" x 13" / 31 cm x 33 cm	1 Jar / Bag 90 / Bag

1 : Although Lila's is completely safe, it is not recommended for use on ophthalmic instruments. Use Lila's Scrub Pad on its own.

2: The stain remover should not be mixed with bleach and is for external use only to remove instrument stains.