

nasal ventilation

Single ventilation tube

Single silicone nasal tubes designed to reduce the stress related to ventilation and to improve the comfort of premature babies and newborns.

- Open tip
- 15mm adaptor



Single ventilation tube

Double tube in silicone for nasal neonatal ventilation

Length of tube on 2 branches: 6 cm, except code 2596.30: 7 cm.

Open and rounded distal tips and centimetric markings.

The proximal end has an ISO standard 15 mm adaptor. This tube is specially designed to avoid tracheal intubation which is very traumatic, and to improve the comfort of premature babies and neonates.

4 sizes are available: Ext. Ø 2.7 mm - 3.4 mm - 4.0 mm - 5.0 mm



Double tube in silicone for nasal neonatal ventilation

Product Code	Int. ø	Ext. ø	Diameter	Length	Units/Case
Single Tube					
2595.15	1.5 mm	2.7 mm	8 Fr	9 cm	50
2595.20	2.0 mm	3.4 mm	10 Fr	9 cm	50
2595.25	2.5 mm	4.0 mm	12 Fr	9 cm	50
2595.30	3.0 mm	5.0 mm	14 Fr	9 cm	50
Double Tube					
2596.15	1.5 mm	2.7 mm	8 Fr	6 cm	50
2596.20	2.0 mm	3.4 mm	10 Fr	6 cm	50
2596.25	2.5 mm	4.0 mm	12 Fr	6 cm	50
2596.30	3.0 mm	5.0 mm	15 Fr	7 cm	50