


naso-gastric

Feeding Tubes

Feeding tubes in DEHP-free PVC (for short-term use) and polyurethane (for medium-term use) with flexible female mount and cap.

- Smoothly rounded, closed distal tip and opposed lateral eyes
- Radiopaque
- Available in various lengths: Centimeter graduations (5cm to 25cm)



Standard tube

PVC tubes

Size	04 Fr	05 Fr	06 Fr	07 Fr	08 Fr	09 Fr	10 Fr	12 Fr	Units/Case
40 cm	310.04	310.05	310.06	310.07	-	-	-	-	50
50 cm	-	-	-	-	310.08	310.09	310.10	310.12	50
75 cm	311.04	311.05	311.06	-	311.08	-	311.10	-	50
125 cm	312.04	-	312.06	312.07	312.08	312.09	312.10	312.12	25

PUR tubes

Size	04 Fr	05 Fr	06 Fr	08 Fr	Units/Case
40 cm	I310.04	I310.05	I310.06	-	25
50 cm	-	-	-	I310.08	25
75 cm	-	I311.05	I311.06	I311.08	25