

# NIOSH APPROVED N95 MODEL L188 PARTICULATE RESPIRATOR

Designed for **high filtration capacity** to protect wearers from airborne particles, this respirator is **hypoallergenic** and **latex free**. The flat folded-style adapts well to wearer's facial features and padded nose piece provides comfort when worn over a prolonged period of time. This N95 respirator **has been approved by NIOSH** (National Institute for Occupational Safety and Health), the global standard and your only choice for protection you can truly trust.



**Performance Standard:** The respirator is Certified under 42 CFR 84 and NIOSH approval "TC-84A-6973"

**Filtration Efficiency:**  $\geq 95\%$

**Product Type:** Disposable Respirator Mask

**Certifications:** Filtration efficiency  $\geq 95\%$ . Inhalation resistance, exhalation resistance, and Sodium Chloride (NaCl) N95 tested and conform to NIOSH criterion. (NIOSH procedure TEB-APR-STP-0007, RCT-APR-STP-003, and TEB-APR-STP-0059). Health Canada compliant. Approved by Health Canada as an authorized medical device for uses related to COVID-19.

## FEATURES

- Conforms highest N95 safety standards, NIOSH N95 (National Institute for Occupational Safety and Health)
- Comfortable fit: flat folded-style adapts well to wearer's facial features and padded nose piece provides comfort when worn over a prolonged period of time

## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
SD-L188-R0	N95 Respirator Mask - Duckbill Style	20 / Box