

NIOSH APPROVED N95 MODEL L288 PARTICULATE RESPIRATOR

Designed for **high filtration capacity** to protect wearers from airborne particles, this N95 respirator is **hypoallergenic, latex free**, and features an adjustable **padded nose piece** for a truly comfortable fit. This N95 respirator **has been approved by NIOSH** (National Institute for Occupational Safety and Health), the global standard and your only choice for protection you can truly trust.



Performance Standard: The respirator is Certified under 42 CFR 84 and NIOSH approval "TC-84A-7228"

Filtration Efficiency: $\geq 95\%$

Certifications: Filtration efficiency $\geq 95\%$. Inhalation resistance, exhalation resistance, and Sodium Chloride (NACL) N95 tested and conform to NIOSH criterion. (NIOSH procedure TEB-APR-STP-0007, RCT-APR-STP-003, and TEB-APR-STP-0059). Health Canada compliant. Approved by Health Canada as an authorized medical device for uses related to COVID-19.

FEATURES

- Conforms to highest N95 safety standards, NIOSH N95 (National Institute for Occupational Safety and Health)
- Comfortable fit: adjustable padded nose piece for comfort when worn over a prolonged period of time

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
SD-L288-R0	N-95 Respirator Mask - Cone Style	20 / Box