


# nutri<sup>line</sup><sup>®</sup>

TWINFLO

## A 2Fr dual lumen catheter with peelable needle

nutri<sup>line</sup><sup>®</sup> twinflo has been designed specifically for babies who require multiple infusions. Its dual lumen tip design offers enhanced fluid management and avoids incompatible drugs mixing.



Product Code	Catheter			Introducer		Stylet	Units / Case
	Length	Priming Case Volume	Flow rate (pressure 1bar)	Type	Size		
I252.232G	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	No introducer	-	No	10
I252.232M	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	Microflash	20 G	No	10
I252.230G	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	Microflash	20 G	Yes	10

## Kits & Accessories

Product Code	Description	Units/Box
7371.19	20 G Microflash Introducer	10
I262.30N	20 G Splitting Needle	20
AMS-850TC	V-Cutter Catheter Trimming Tool	20



Microflash



Splitting needle



V-Cutter