

nutriSafe2®

A complete range dedicated to safe enteral nutrition. Designed to be fail-safe, the nutriSafe2® system incorporates a unique connection that locks securely so it won't slip apart and yet cannot be connected to a standard 6% tapered luer fitting.

Safe - The risk of misconnection errors is significantly reduced for patient safety.

Secure - Accidental disconnection is nearly impossible, ensuring precious nutrients are received.

Simple - It does not change the technique of the users and therefore does not require any special training.

nutriSafe2® Feeding tubes

nutriSafe2® feeding tubes are available in a range of sizes, lengths and materials to suit patient needs. Numbered graduations every centimeter aid placement and monitoring, and the rounded tip provides atraumatic insertion.

PVC (DEHP-free): for use up to 5 days.

Polyurethane (PUR): for use up to 30 days.

Silicone: for use up to 30 days.

The nutriSafe2® connection allows the use of these tubes only with other nutriSafe2® devices..

PVC tubes



Size	04 Fr ●	05 Fr ●	06 Fr ●	08 Fr ●	10 Fr ●	Units/Case
40 cm	364.042	364.052	364.062	364.082	-	25
50 cm	361.042	361.052	361.062	361.082	361.102	50
75 cm	363.042	363.052	363.062	363.082	363.102	25
125 cm	362.042	362.052	362.062	362.082	362.102	25

PUR tubes



Size	04 Fr ●	05 Fr ●	06 Fr ●	08 Fr ●	10 Fr ●	Units/Case
40 cm	1364.042	1364.052	1364.062	1364.082	-	25
50 cm	1361.042	1361.052	1361.062	1361.082	1361.102	50
75 cm	1363.042	1363.052	1363.062	1363.082	1363.102	25
125 cm	1362.042	1362.052	1362.062	1362.082	1362.102	25

Silicone tubes



Size	04 Fr ●	05 Fr ●	06 Fr ●	08 Fr ●	10 Fr ●	Units/Case
50 cm	2331.042G	2331.052G	2331.062G	2331.082G	2331.102	10
125 cm	2332.042G	-	2332.062G	2332.082G	2332.102	10
125 cm (weighted)	-	-	2395.062G	2395.082G	-	10

nutri^{safe}2[®] Syringes

A complete range of syringes, from 1 ml to 60ml for administration of enteral nutrition and oral medication.

Amber colored syringes allow preservation of light sensitive drugs.
They are available with caps.

These syringes only connect with other items in the nutri^{safe}2[®] range.



Transparent syringes



Amber syringes

C-GON - Transparent syringes

Product Code	Volume	Graduations	Connection	Units/Case
I015.012	1 ml	0.01 ml	Centered / slip	100
I015.022	2.5 ml	0.1 ml	Centered / slip	100
I015.052	5 ml	0.2 ml	Centered / slip	100
I015.102	10 ml	0.2 ml	Centered / lock	100
I015.202	20 ml	1 ml	Centered / lock	50
I015.352	35 ml	1 ml	Centered / lock	50
I015.602	60 ml	1 ml	Offset / lock	25

Monoject - Transparent syringes

Product Code	Volume	Graduations	Connection	Units/Case
I015.062M	6 ml	0.1 ml	Centered / slip	100
I015.202M	20 ml	1 ml	Centered / lock	50
I015.352M	35 ml	1 ml	Offset / lock	25
I015.602M	60 ml	1 ml	Offset / lock	25

Amber syringes (bulk non-sterile)

Product Code	Volume	Graduations	Connection	Units/Case
I022.012	1 ml	0.01 ml	Centered / slip	50
I022.022	2.5 ml	0.1 ml	Centered / slip	50
I022.052	5 ml	0.2 ml	Centered / lock	50
I022.102	10 ml	1 ml	Centered / lock	50

nutrisafe²® Extension lines

Extension lines to allow for positioning of the syringe driver away from the incubator and to move handling away from the newborn.

Available in standard bore for breast milk or in large bore for fortified milk.



Product Code	Int. ø	Length	Clamp	Units/Case
368.032-50	1.5 mm	30 cm	No	50
368.152-50	1.5 mm	150 cm	No	50
369.152	2.8 mm	150 cm	No	25

nutrisafe²® Caps

A range of caps to secure the contents of vials containing oral / enteral medications. These caps allow access to the medication only from a nutrisafe²® syringe, to avoid risk of accidental injection through IV systems.

Caps for syringes and tubes prevent contamination of the feed or drugs between feeding sessions.

Needleless vial access caps

Product Code	Connector	Int. ø	Units/Case
821.022	male	20 mm	25



Needleless vial access caps

Caps for syringes and tubes

Product Code	Description	Units/Case
828.01	Cap for syringe (x1)	100
828.08	Cap for syringe (x8)	100



Syringe cap

Vented caps for syringes - for gravity feeding

Product Code	Description	Units/Case
818.22	20 ml	50
818.62	60 ml	50



Vented cap

nutri^{safe}2[®] sampling straws



nutri^{safe}2[®] drawing-up straws allow the filling of nutri^{safe}2[®] syringes with feed or medication from all sizes of containers and bottles.

Product Code	Size	Bore	Units/Case
817.052	5 cm	Standard bore	50
817.152	15 cm	Standard bore	50
817.153	15 cm	Large bore	50
817.202	20 cm	Standard bore	50
817.203	20 cm	Large bore	50
817.302	30 cm	Large bore	50
817.502	50 cm	Standard bore	20

nutri^{safe}2[®] accessories

Connectors, stopcocks and accessories to facilitate daily use of the nutri^{safe}2[®] range.

Connectors

Product Code	Description	Units/Case
368.02	Adaptor for feeding sets	100
368.12	Taper connector	30
368.32	Extension tube with connector for suction	25
368.42R	Y connector with oral port	100



Adaptor for feeding sets



Y connector



Taper Connector



Extension Tube with connector

Stopcock

Product Code	Description	Units/Case
5802.01	3-way stopcock	50



Stopcock

Adapters

Product Code	Description	Units/Case
828.88	Oral tip adapter	1000
828.99	Medicine cap adapter	1000



Oral tip adapter



Medicine cap adapter