


# ASEPT®

## Pleural & Peritoneal Drainage System

### ASEPT Catheter and Insertion Kits

The ASEPT Pleural and Peritoneal Drainage Systems are designed for long-term, intermittent drainage of symptomatic, recurrent fluid accumulation.



- |  |                                     |                         |                       |
|--|-------------------------------------|-------------------------|-----------------------|
| 1. ASEPT Catheter  | 6. 16G x 1.0" (25.4mm) Needle       | 12. Introducer Needle   | 18. Tunneler          |
| 2. ASEPT Drainage Line   | 7. 25G x 1.0" (25.4mm) Needle       | 13. Tearaway Introducer | 19. Forceps           |
| 3. 5-in-1 Drain Line Adapter                                       | 8. 22G x 1.5" (38.1mm) Needle       | 14. 8F (2.67mm) Dilator | 20. Safety Scalpel    |
| 4. Chloraprep® 10.5mL Vial   | 9. Tegaderm® Self Adhesive Dressing | 15. 12F (4.0mm) Dilator | 21. Fenestrated Drape |
| 5. Synthetic Suture, 2-0 Straight and Synthetic Suture, 3-0 Curved | 10. Gauze Pads, 4"x4"               | 16. 10mL Syringe        |                       |
|  | 11. Foam Catheter Pad               | 17. J-Tip Guidewire     |                       |


The primary component of the ASEPT Drainage system is the catheter, which has an aseptic, self-sealing valve attached that provides a secure and reliable connection. The valve can be connected to an ASEPT vacuum drainage bottle, drainage bag, or a stand-alone drainage line to connect to other fluid collection equipment.

#### Features:

- Secure connection
- Cap-less and Aseptic design
- Needle-free and swabable
- Easy to replace
- Self-sealing

### ASEPT Drainage Bottles and Bags

The ASEPT Drainage kit configuration includes a 600mL or 1000mL vacuum bottle set, a 2000mL gravity drainage bag, and the procedure pack with all of the dressing change items needed to complete an at-home drainage procedure. A standalone drain line is available to connect to the ASEPT drainage catheter for drainage with other collection equipment during in-patient procedures.



#### Bottle Features:

- Convenient vacuum release clamp
- Green vacuum indicator
- Luer-lock connector



#### Bag Features:

- Large 2000mL volume
- Fluid release valve for ease of bag emptying

## ASEPT Accessories

PFM Medical, Inc. is proud to offer a full line of accessories ensuring that the physician or patient will have everything or she needs to complete the drainage procedure. Our accessory portfolio contains the following items:

### ASEPT Replacement Valve



### ASEPT Drainage Line Set



### ASEPT Procedure Pack



## Ordering Information

### ASEPT Pleural/Peritoneal Insertion Kits

| REF       | Description                       | Box Qty (ea) |
|-----------|-----------------------------------|--------------|
| M7001-CAN | Pleural Drainage Insertion Kit    | 1            |
| M7008-CAN | Peritoneal Drainage Insertion Kit | 1            |

### ASEPT Drainage Bottles & Bags

|       |  |    |
|-------|--|----|
| M7002 | 600 mL Drainage Kit*                                     | 10 |
| M7080 | 600 mL Drainage Starter Kit**                            | 5  |
| M7072 | 600 mL Bottle Only                                       | 10 |
| M7000 | 1000 mL Drainage Kit*                                    | 10 |
| M7081 | 1000 mL Drainage Starter Kit**                           | 5  |
| M7070 | 1000 mL Bottle Only                                      | 10 |
| M7060 | 2000 mL Drainage Bag Starter Kit** (Peritoneal Use Only) | 10 |
| M7058 | 2000 mL Drainage Bag (Peritoneal Use Only)               | 10 |

### ASEPT Accessories

|       |  |    |
|-------|--|----|
| M7040 | Procedure Pack (Traditional Dressing Change) | 10 |
| M7004 | Catheter Replacement Valve                   | 5  |
| M7005 | Drainage Line Set (includes 5-in-1 adapter)  | 10 |

\*ASEPT Drainage Kit includes bottle and traditional dressing change procedure pack.

\*\* ASEPT Drainage Starter Kit recommended for patient discharge. Includes bottles & traditional dressing change procedure packs.

## Indications for Use

### ASEPT Pleural Indications

The ASEPT Pleural Drainage System is intended for long term, intermittent drainage of symptomatic, recurrent, pleural effusions, including malignant pleural effusions and other pleural effusions that do not respond to treatment of the underlying disease.

### ASEPT Peritoneal Indications

The ASEPT Peritoneal Drainage System is indicated for periodic drainage of recurrent and symptomatic malignant ascites. The catheter is intended for long term access of the peritoneal cavity in order to relieve symptoms such as dyspnea.