


SMOOTH MOVER



PRODUCT DESCRIPTION

An economical ultra-smooth polyethylene patented board especially designed to furnish expert patient transference with minimal trauma to the patient. The Smooth Mover is capable of safely moving extremely large patients with minimal staffing requirements. It is radiotranslucent; it can be left in place during x-ray procedures. Is lightweight yet durable. Perfect for daily patient moving needs: EMS or hospital, nursing home or clinic applications. Easy on ...easy off. Yes, even one person can move a large patient!

Optional Accessory:

Pull Straps... loops straps slip through hand-holes enabling easy board pick-up on hard surfaces tables (x-ray, etc.) and minimize reach required when transferring patients from bed to stretcher to table during routine moving.

PRODUCT ORDERING INFORMATION

Hospital Shifter Model SM3

Is widely used for patient shifting, bed to stretcher to table, etc..
Dimensions: 3/16" x 22" x 72"
Weight: 9 lbs.

Ambulance Patient Shifter Model SM40

Is used as an aid in moving patients from vehicle to stretcher or backboard.
Dimensions: 3/16" x 18" x 48"
Weight: 5½ lbs.