



TABLES, CABINETS, CARTS AND STANDS

DELIVERING TOTAL HEALTHCARE VALUE SINCE 1935



CLOSED CASE CARTS

Our Closed Case Carts are built to last and come in a wide range of sizes with flexible shelving options, ergonomic door assemblies, smart design choices, important safety features, superior maneuverability, and optional accessories. They are highly robust and durable yet lightweight thanks to our one-piece body-wrap design of 18-gauge superior stainless steel.

All closed case carts are backed by an outstanding **7 year warranty** on structural and all components (excluding casters) with one year labor included.

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS (L x W x H) IN / MM	CAPACITY (CU. FT.) *
OR-3810	SINGLE DOOR / STANDARD HEIGHT	Shelf: 24" x 24" (610 mm x 610 mm) Cabinet: 27" x 25 1/2" x 33" (686 mm x 648 mm x 838 mm) Overall: 33 1/4" x 28" x 39 1/2" (845 mm x 711 mm x 1003 mm)	11.5
OR-3815	SINGLE DOOR / TALL DESIGN	Shelf: 24" x 24" (610 mm x 610 mm) Cabinet: 27" x 25 1/2" x 48" (686 mm x 648 mm x 1219 mm) Overall: 33 1/4" x 28" x 54 1/2" (845 mm x 711 mm x 1384 mm)	16.5
OR-3820	DOUBLE DOOR / STANDARD HEIGHT	Shelf: 36" x 24" (914 mm x 610 mm) Cabinet: 39" x 25 1/2" x 33" (991 mm x 648 mm x 838 mm) Overall: 45 1/4" x 28" x 39 1/2" (1149 mm x 711 mm x 1003 mm)	16.6
OR-3825	DOUBLE DOOR / TALL DESIGN	Shelf: 36" x 24" (914 mm x 610 mm) Cabinet: 39" x 25 1/2" x 48" (991 mm x 648 mm x 1219 mm) Overall: 45 1/4" x 28" x 54 1/2" (1149 mm x 711 mm x 1384 mm)	24.2

ACCESSORIES

ITEM #	DESCRIPTION
OR-3875	CARD HOLDER 6" W x 4" H (152mm x 102mm)
OR-3876	CARD HOLDER 9" W x 11" H (229mm x 279mm)
OR-3877	SECURITY / RING-SEAL BRACKETS
OR-3881	PUSH / PULL HANDLE
OR-3871-S	PADDLE LATCH ASSEMBLY SMALL
OR-3871-L	PADDLE LATCH ASSEMBLY LARGE
OR-3840	SOLID SLIDE-OUT SHELF 24" W x 24" H - Only for single door models
OR-3841	PERFORATED SLIDE-OUT SHELF 24" W x 24" H - Only for single door models
OR-3843	WIRE SLIDE-OUT SHELF 24" W x 24" H - Only for single door models
OR-3845	SOLID SLIDE-OUT SHELF 24" W x 36" H - Only for double door models
OR-3846	PERFORATED SLIDE-OUT SHELF 24" W x 36" H - Only for single door models
OR-3848	WIRE SLIDE-OUT SHELF 24" W x 36" H - Only for single door models



OR-3825

KEY FEATURES

FLEXIBLE SHELVING

- Adjustable open shelf design provides additional shelf width and depth for large containers
- Shelves types include solid, perforated, or wire with optional roll-out featuring integral stops to prevent tipping during use

ERGONOMIC DOOR ASSEMBLIES

- Strong double pan design features bottom edge drainage holes (outside positioning)
- Heavy gauge hinge assemblies are easily removed for cleaning/adjustment
- SteriLock spring loaded door latch prevents accidental opening during transportation and ensures a fully open position during washing, loading, and unloading



SAFETY FEATURES

- Full perimeter bumper made from high temperature resistant rubber provides impact protection
- Permanently sealed seams and crevices prevents water accumulation and sources of contamination

SMART DESIGN

- Top surface features an integral drip-proof ledge to ensure proper drainage of accidental spills without contaminating cart interior
- Cabinet bottom pan is sloped 5 degrees and chassis support frame is perforated to facilitate drainage
- Externally mounted integrally welded push-pull handle maximizes interior storage space

SUPERIOR MANEUVERABILITY

- Each cart features two all-swivel, one steer-lock, and one total-lock caster to facilitate easy steering, handling, and positioning
- Stainless steel caster body is cart washer compatible and provides a high load rating of 225 lbs per wheel
- Optional retractable linkage/training assembly provides for easy manual cart training/towing

OPTIONAL ACCESSORIES

- Choose amongst optional card holders, paddle latch assembly, linkage, and roll-out shelving



OPEN CASE CARTS

Our Open Case Carts are built to last using Series 300 stainless steel with a #4 polished finish and come in two versatile sizes and three heights with flexible shelving options, smart design choices, important safety features, and superior maneuverability. The tubular design offers maximum rigidity, and the frame assembly is structurally supported by a heavy gauge fully welded chassis assembly for maximum load bearing capacity.

All open case carts are backed by an outstanding **7 year warranty** on structural and all components (excluding casters) with one year labor included.

KEY FEATURES

FLEXIBLE SHELVING

- Adjustable open shelf design provides additional shelf width and depth for large containers and can be easily removed/repositioned on the cross supports in 5" (127 mm) increments
- Shelves are available in solid or perforated stainless steel and are die-formed for maximum load bearing capacity

SMART DESIGN

- Work height carts have a fully welded solid top surface with perimeter channel supports for structural integrity
- Externally mounted integrally welded push-pull handle maximizes interior storage space
- All welds and joints are polished and blended for a smooth finish

SAFETY FEATURES

- Full perimeter bumper made from high temperature resistant rubber provides impact protection
- Standard corner bumpers are available as an option
- Permanently sealed seams and crevices prevents water accumulation and sources of contamination
- All surfaces are smooth and free from sharp edges that may result in sterile pack damage or worker injury



SUPERIOR MANEUVERABILITY

- Each cart features two all-swivel, one steer-lock, and one total-lock caster to facilitate easy steering, handling, and positioning
- Stainless steel caster body is cart washer compatible and provides a high load rating of 225 lbs per wheel
- Optional retractable linkage/training assembly provides for easy manual cart training/towing





PRODUCT SELECTION

FEATURES	OR-3920 & OR-3921-PS (add suffix A, B or C to item # to adjust features' specifications)		
	-A	-B	-C
OVERALL HEIGHT	39 1/2" (1003 mm)	39 1/2" (1003 mm)	54" (1372 mm)
CAPACITY (CU. FT)	16.6	23.9	28.9
SHELF SUPPORT POSITION	3	3	3
CASTERS SIZE	5" (127 mm)	5" (127 mm)	8" (203 mm)
OVERALL DIMENSION	45 1/4" L x 28" W (1149 mm x 711 mm)		
SHELF DIMENSION	36 1/4" L x 24" W (921 mm x 610 mm)		

Item # OR-3921-PS has perforated shelves

ACCESSORIES

ITEM #	DESCRIPTION
OR3982	BUMPER ASSEMBLY FLUSH MOUNT



UTILITY CARTS

Our Utility Carts are built to last using heavy-gauge grade 18-8 stainless steel and come in four versatile models for general purpose transport/storage use of supplies, instruments, or equipment in the Operating Room, trauma, emergency or minor procedure area.

All utility carts are backed by our **1 year warranty** on structural and all components (excluding casters).

KEY FEATURES

SMART DESIGN

- Fully welded and polished to a #4 finish
- Spill-proof 1 inch (51 mm) lip on three sides of top surface and middle shelf
- Non-marking bumpers on push/pull handle and front corners of the cart
- **Reinforcing steel rod** is secured below the top shelf for **added strength and stability**
- **Steel plate** is affixed underneath all shelves to prevent wobble and bottoming-out

SUPERIOR MANEUVERABILITY

- Each cart features high load rating swivel casters that are mounted to heavy-duty chassis base frame
- Push/pull handle runs the full width of the cart and is positioned for ease of use



OR-3692F

PRODUCT SELECTION

ITEM #	SHELF SIZE (W x L)	SHELF CLEARANCE	OVERALL HEIGHT	CASTER DIAMETER	LOAD CAPACITY
OR-3692A	15 1/2" x 24" (394 mm x 610 mm)	12" (305 mm)	32" (813 mm)	3" (76 mm)	200 lbs (91 kg)
OR-3692B	17 3/4" x 30" (451 mm x 762 mm)	11 3/4" (298 mm)	32" (813 mm)	3" (76 mm)	200 lbs (91 kg)
OR-3692D	17 3/4" x 30" (451 mm x 762 mm)	12 1/8" (308 mm)	33 1/4" (845 mm)	4" (102 mm)	400 lbs (182 kg)
OR-3692F	20 3/4" x 35" (527 mm x 889 mm)	13 1/2" (343 mm)	37 1/8" (943 mm)	5" (127 mm)	500 lbs (224 kg)
OR-3692-WO*	19" x 31" (483 mm x 787 mm) Area 4.2 SQ FT (0.4 SQ M) Volume: 12 CU FT (0.3 CU M) Weight: 42 lbs. (19.1 kg)	Customizable	34" (864 mm)	4" (102 mm)	500 lbs (181 kg)

*Three sides have 1" high raised edge on top and middle shelf.

NOTE: Due to continual product development, specifications are subject to change without notice.

ANESTHESIA CARTS

Our Anesthesia Cart is a flexible storage system for the various anesthesia supplies, medications, and instrumentation required in the operating room. Built to last in stainless steel 304 with a #4 polish, it is easy to clean and maintain.

All anesthesia carts are backed by our **1 year warranty** on structural and all components (excluding casters).

KEY FEATURES

SMART DESIGN

- Comes complete with six full extension drawer assemblies in various depths: 3 x 3" (76 mm), 2 x 6" (152 mm), 1 x 9" (228 mm)
- Smooth surface with a three-sided guard rail
- Over bridge assembly with optional accessory bins and/or shelf provides easily accessible additional storage space

SUPERIOR MANEUVERABILITY

- Caster assemblies with high load bearing capacity and grey non-marking rubber tires provide superior maneuverability
- Full perimeter bumper assembly protects cart body during transport

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS (L x W x H)
OR-3680	ANAESTHESIA CART	27" X 22" X 40" (686 mm x 559 mm x 1016 mm)



OR-3680

Products shown in Anaesthesia Cart not included

ACCESSORIES

ITEM #	DESCRIPTION
OR-3680-DV3	DIVIDER KIT for 3" / 76 mm drawer
OR-3680-DV6	DIVIDER KIT for 6" / 152 mm drawer
OR-3680-GH	GLOVE BOX HOLDER
OR-3680-NC3	NARCOTIC CABINET WITH LOCK inside 3" / 76 mm drawer
OR-3680-NC6	NARCOTIC CABINET WITH LOCK inside 6" / 152 mm drawer
OR-3680-OR	OVERHEAD RAIL ASSEMBLY FOR BINS Complete with 2 lateral supports
OR-3680-PB	POWER BAR, HOSPITAL GRADE Complete with 10 ft / 3 M cord
OR-3680-SH	OVERHEAD SHELF ASSEMBLY 26" L x 10" W / 660 mm L x 254 mm W
OR-3680-SR	SIDE RAIL BRACKET
OR-3680-VS	VELCRO STRAP / BRACKET KIT



INSTRUMENT/ SUPPLY CABINET

Our Instrument/Supply Cabinets are designed to provide enclosed storage for O.R. supplies, medications, instrumentation, and apparatus.

All Instrument/Supply Cabinets are backed by our **1 year warranty** on structural and all components.

KEY FEATURES

- Cabinets come with four easily adjustable shelves (on 1" / 25 mm increments) and are fully interchangeable with other instrument/supply cabinets of similar width
- Standard sliding glass doors are framed in stainless steel and provide quick access to cabinet contents (optional door systems available)
- Cabinets are 22" / 559 mm deep x 84" / 2134 mm high at front
- Optional sloped top section adds 9 ½" / 241 mm to overall cabinet height
- Optional depth cabinets of 17 ¾" / 450 mm are available upon request



OR-3782

PRODUCT SELECTION

ITEM #	DESCRIPTION
OR-3782-A	INSTRUMENT/SUPPLY CABINET Width: 47" (1194 mm)
OR-3782-B	INSTRUMENT/SUPPLY CABINET Width: 35" (889 mm)
OR-3782-C	INSTRUMENT/SUPPLY CABINET Width: 24" (610 mm)
OR-3782-GH	HINGED GLASS DOORS
OR-3782-RF	RECESSED FLANGE KIT
OR-3782-ST	SLOPED TOP SECTION

INSTRUMENT TABLES

Our Instrument Tables come in six versatile models and are an economical back table for instruments or dressings. Fully welded heavy gauge stainless steel construction polished to a #4 finish provides excellent rigidity and extended product life.

All instrument tables are backed by an outstanding **7 year warranty** on structural and all components (excluding casters).

KEY FEATURES

- Channel support sub-frame for maximum weight capacity and structural rigidity
- Fully welded seamless tubular legs H-frame supports or lower shelf assembly
- Non-marking swivel design casters come in 2" (51 mm) or 3" (76 mm) depending on table size
- All models are 34" (864 mm) high and feature a lower shelf
- Open tubular base design is available upon request (reference # OR-3643)



PRODUCT SELECTION

ITEM #	DIMENSIONS (L x W)
OR-3641A	20" x 16" (506 mm x 406 mm)
OR-3641B	30" x 16" (762 mm x 406 mm)
OR-3641C	34" x 18" (864 mm x 457 mm)
OR-3641D	36" x 20" (914 mm x 508 mm)
OR-3641H*	46" x 24" (1168 mm x 610 mm)
OR-3641L*	58" x 24" (1473 mm x 610 mm)

*H and L models supplied with 3" / 76 mm caster assemblies

ACCESSORIES

ITEM #	DESCRIPTION
OR-3641X-R3	THREE-SIDED GUIDE RAIL Assembly for Table Top or Lower Shelf (Specify size)
OR-3641X-R4	FOUR-SIDED GUIDE RAIL Assembly for Table Top or Lower Shelf (Specify size)
OR-3640	DRAWER ASSEMBLY Inside Dimensions: 15" L x 16 3/4" W x 3 3/4" D (381 mm x 426 mm x 95 mm)

EQUIPMENT TABLES

Our Equipment Tables come in three sizes and are a sterile and efficient solution for organizing and storing equipment and supplies. Fully welded heavy gauge stainless steel construction is easy to clean and maintain.

All equipment tables are backed by our **1 year warranty** on structural and all components (excluding casters).

KEY FEATURES

- All models include stainless steel welded drawer liners, a lower shelf, and three -sided ribbon rail
- Heavy-duty conductive swivel casters come in 2" (51 mm) or 3" (76 mm) depending on table size



PRODUCT SELECTION

ITEM #	DESCRIPTION (L x W x H)
OR-3644	EQUIPMENT/INSTRUMENT TABLE WITH TWO DRAWERS, GUARD RAIL, AND LOWER SHELF 36" x 18" x 34" (914 mm x 457 mm x 864 mm)
OR-3648	EQUIPMENT/INSTRUMENT TABLE WITH ONE DRAWER, GUARD RAIL, AND LOWER SHELF 20" x 18" x 34" (508 mm x 457 mm x 864 mm)
OR-3650	ANESTHETIST'S TABLE WITH THREE DRAWERS (2 X 4"/102MM DEEP, 1 X 8"/204MM DEEP), GUARD RAIL, AND LOWER SHELF 20" x 18" x 34" (508 mm x 457 mm x 864 mm)

MAYO STANDS

Our Mayo Stands come in seven models and are a great choice for general clinical use and operating room applications. The stainless steel seamless tubular construction is fully welded and polished for ease of cleaning and maximum structural support.

All mayo stands are backed by an outstanding **7 year warranty** on structural and all components (excluding casters).


KEY FEATURES

- Auto-clutch height adjustment system automatically locks tray at desired height and provides for smooth raising and lowering
- Safety mechanism prevents the tray frame from unexpectedly dropping
- Casters are ball bearing swivel 2" (51 mm) diameter
- All models feature a tray height adjustment range of 38" to 58" (965 mm to 1473 mm)

PRODUCT SELECTION

ITEM #	DESCRIPTION	TRAY SIZE (L X W)
OR-3657	EMERGENCY MAYO STAND WITH 4 CASTER BASE, AND T-LOCK Manually Operated	19"x13" (483 mm x 330 mm)
OR-3658	STANDARD TRAY Push-Button Sled-Type Base and Auto-Clutch System	19"x13" (483 mm x 330 mm)
OR-3662	LARGE TRAY Push-Button Sled-Type Base and Auto-Clutch System	21"x16" (533 mm x 406 mm)
OR-3660	STANDARD TRAY Manually Operated with Two Casters Sled-Type Base and T-Lock Handle	19"x13" (483 mm x 330 mm)
OR-3661	LARGE TRAY Manually Operated with Two Casters Sled-Type Base and T-Lock Handle	21"x16" (533 mm x 406 mm)
OR-3665	STANDARD TRAY Foot-Operated Mayo Stand Sled-Type Base and Foot Release Auto-Clutch System	19"x13" (483 mm x 330 mm)
OR-3666	LARGE TRAY Foot-Operated Mayo Stand Sled-Type Base and Foot Release Auto-Clutch System	21"x16" (533 mm x 406 mm)

STAINLESS STEEL TRAYS ARE POLISHED, AND DIE-FORMED WITH RADIUS CORNERS. THEY COME IN TWO STANDARD SIZES



OR-3665

STRADDLE (NEUROSURGERY) TABLES

Our Straddle (Neurosurgery) Tables come in four lengths and are a convenient height-adjustable instrument table for cases involving many instruments and/or diagnostic equipment. Heavy-duty polished stainless steel construction with fully welded tubular support base.

All straddle (neurosurgery) tables are backed by our **1 year warranty** on structural and all components (excluding casters).

KEY FEATURES

- Direct-action drive mechanism provides smooth adjustment with no creeping or slippage
- Optional accessories include stainless steel extension shelves and marine-edge spill proof stainless steel work surface
- Feature four 3" (76 mm) grey non-marking casters (2 with swivel, 2 with brake)
- Height-adjustable from 36" to 58" (914 mm to 1473 mm) with stainless steel collapsible crank handle
- All models feature a standard 27.5" (700 mm) minimum clearance between the vertical leg supports



OR-3646X

Standard 27 1/2" / 700 mm minimum clearance between the vertical leg supports

PRODUCT SELECTION

ITEM #	DESCRIPTION
OR-3646X	STRADDLE TABLE, WITH STANDARD POSITION CRANK HANDLE Available lengths: A: 24" / 610 mm B: 32" / 813 mm C: 40" / 1016 mm D: 46" / 1168 mm

Add suffix to Model # for work surface length required.
Standard Table / Extension Shelf width is 34" / 864 mm.

Stainless steel extension shelves and marine-edge work surfaces are available as optional features.

ACCESSORIES

ITEM #	DESCRIPTION
OR-3646ES	STAINLESS STEEL EXTENSION SHELF Available lengths: 4" and 8" / 102 mm and 203 mm depending on table width
OR-3646LH	LOWER POSITION CRANK HANDLE ASSEMBLY
OR-3646ME-X	MARINE-EDGE SPILL PROOF STAINLESS STEEL WORK SURFACE


LINEN INSPECTION TABLES

Our Linen Inspection Tables provide a large, LED illuminated work surface for the inspection of linen drapes and packing wraps. Available in 2 sizes and constructed from stainless steel with a sealed translucent plexiglass work surface. The stainless steel frame and legs ensures durability for extended product life.

All linen inspection tables are backed by our **1 year warranty** on structural and all components (excluding casters).

KEY FEATURES

- LED fixtures provide bright, consistent lighting, are power efficient and last longer than traditional bulbs
- Optional lower shelf provides extra storage space
- Tubular legs are available with optional heavy-duty grey non-marking casters with brakes
- Features height adjustable levelers
- CSA Approved hospital grade power cord (10 feet / 3 meters) and light switch
- Easily accessible lamp assembly for bulb replacement
- Service required: 120V, 60Hz



SP-8210 LS

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS (L x W x H)
SP-8210 LS	LINEN INSPECTION TABLE WITH LOWER SHELF (Standard Length)	72"x 32"x 36" (1829 mm x 813 mm x 914 mm)
SP-8210	LINEN INSPECTION TABLE (Standard Length)	72"x 32"x 36" (1829 mm x 813 mm x 914 mm)
SP-8215	CASTER KIT FOR STANDARD TABLE (Set of Four) Factory Installed	
SP-8220	LINEN INSPECTION TABLE (Extended Length)	92"x 32"x 36" (2337 mm x 813 mm x 914 mm)
SP-8220 LS	LINEN INSPECTION TABLE WITH LOWER SHELF (Standard Length)	92"x 32"x 36" (2337 mm x 813 mm x 914 mm)
SP-8225	CASTER KIT FOR EXTENDED TABLE (Set of Six)	

UTILITY INSPECTION TABLE

Our Utility Inspection Tables facilitate the assembly, inspection, and packing of instruments and/or equipment before or after the decontamination and sterilization process. All table components are built to last using heavy gauge Series 300 Grade 18-8 stainless steel and are die-formed, fully welded, and polished to a #4 blended finish.

All Utility Inspection Tables are backed by our **1 year warranty** on structural and all components.

KEY FEATURES

- Rigid channel support design and tubular leg cross bracing ensures maximum stability, high load capacity, and extended product life
- Countertop is supported with die-formed stainless steel channels and integrally welded to the tubular leg frame base assembly
- Tables feature height adjustable stainless steel bullet feet
- Optional drawer assemblies, full length lower shelf, and overshelf assemblies provide flexible storage configurations
- Optional drawer assemblies are die-formed stainless steel with a double pan drawer front and heavy-duty ball bearing slide assemblies
- Custom sizes and designs are available upon request



SP-8724

PRODUCT SELECTION

ITEM #	DESCRIPTION
SP-8710	UTILITY INSPECTION TABLE WITH TUBULAR LEG ASSEMBLY Length available: 44" / 1118 mm
SP-8720	UTILITY INSPECTION TABLE WITH TUBULAR LEG ASSEMBLY Length available: 60" / 1524 mm
SP-8730	UTILITY INSPECTION TABLE WITH TUBULAR LEG ASSEMBLY Length available: 76" / 1930 mm
SP-8740	UTILITY INSPECTION TABLE WITH TUBULAR LEG ASSEMBLY Length available: 92" / 2337 mm
SP-8712	OVERSHELF ASSEMBLY WITH MOUNTING KIT Length available: 44" / 1118 mm
SP-8722	OVERSHELF ASSEMBLY WITH MOUNTING KIT Length available: 60" / 1524 mm
SP-8732	OVERSHELF ASSEMBLY WITH MOUNTING KIT Length available: 76" / 1930 mm
SP-8742	OVERSHELF ASSEMBLY WITH MOUNTING KIT Length available: 92" / 2337 mm
SP-8714	MODULAR DRAWER WITH FOUR SIDED APRON ASSEMBLY Length available: 44" / 1118 mm
SP-8724	MODULAR DRAWER WITH FOUR SIDED APRON ASSEMBLY Length available: 60" / 1524 mm
SP-8734	MODULAR DRAWER WITH FOUR SIDED APRON ASSEMBLY Length available: 76" / 1930 mm
SP-8744	MODULAR DRAWER WITH FOUR SIDED APRON ASSEMBLY Length available: 92" / 2337 mm



SURGMED GROUP

2165 46th Ave., Lachine, Quebec H8T 2P1
Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax: +1.514.631.9083
Email: info@surgmed.com | Website: www.imperialsurgical.com