


Tampa Infusion Catheter

For sonohysterography - the simple, cost-effective choice for ultrasound examination of the uterus.

- Unique introducer sheath (33 cm long) facilitates insertion into cervical Os with minimal trauma
- 7 cm marker provides optimal placement of catheter tip in uterus
- Balloonless configuration is ideal for postmenopausal patients with stenotic cervical canals, as well as for many nulliparous patients
- Small catheter size allows easy placement next to ultrasonography probe



| Item # | Description |
|---------|--|
| 61-2005 | Tampa Catheter (20/Box) 5 French, clear, flexible catheter with introducer sheath, 33 cm long, sterile |

