

Bronchitrac® “L” Set

There is no need to twist the catheter or patient's body. Simply advance the Bronchitrac® “L” through a properly placed endotracheal or tracheostomy tube. The dual-plane, patented curvature of the catheter ensures consistent entry into the left main-stem bronchus.

The time savings and reduced patient trauma provided by the Bronchitrac® “L” make it the obvious choice when suctioning the left-bronchial airway.

- Ideal for suctioning of post-cardiac surgical, neuro and critical care patients
- Can be used for deep-lobe suctioning
- Completely latex-free



Item #	Description
41-7814	21" long, 14 Fr Bronchial Suction Catheter, Slide-cover Connector and (2) Vinyl Gloves, Sterile (20/Box)

TRAKE-fit™ Collar Set

Tracheostomy Tube Holder For Maximum Patient Comfort

- Available in single- or double-collar sets
- Installs in seconds, remains secure
- Maximum comfort with cotton-lined neckband
- No adhesive, no narrow tape
- Completely adjustable
- Narrow tabs available for pediatric and small adult use
- Third-party reimbursable



Item #	Description
43-1850	Double Collar Set (25/Box) Each Set includes: (2) neckbands, (2) wide hook/loop tabs
43-1851	Single Collar Set (25/Box) Each Set includes: (1) neckband, (2) wide hook/loop tabs
43-1862	Double Collar/Narrow Tab Set (25/Box) Each Set includes: (2) neckbands, (2) narrow hook/loop tabs

