

# TransWarmer®

## Transport Mattress

Consistent, Safe Warming for Newborns



The TransWarmer® Infant Transport Mattress provides safe and reliable heat for up to two hours when cold stress is a concern. This thermostable, nontoxic gel mattress nests the infant in soothing warmth, while minimizing vibration during transport. Based on a 24 °C/75 °F start temperature, the mattress reaches a peak temperature of 40 °C/104 °F in less than 60 seconds - quickly improving the infant's admission temperature and reducing the risk of hypothermia.



### Features and Benefits

- Comfortable, gel-filled mattress cradles baby in therapeutic warmth
- Gently conforms to the patient's shape for better nesting
- Minimizes the risk of hypothermia in cesarean deliveries
- Better heat distribution increases blood flow
- Disposable heat pack minimizes risk of cross contamination and infection
- Patented disc mechanism is easy to activate
- TransWarmer is nontoxic, latex- and PVC-free
- Use in labour and delivery, NICU, nursery, OR, ER, and transport units
- 10" wide x 16" tall

Item #	Description
20421	TransWarmer® Infant Transport Mattress (6/Box)