

Gastric Lavage

# TUM-E-VAC® GASTRIC LAVAGE KIT

LATEX FREE | SINGLE USE | DISPOSABLE

Ethox



- Pre-assembled kits are quick and easy to use, reducing the time for lavage procedures
- Promotes quick, efficient lavage of stomach contents
- Closed system helps prevent exposure to stomach contents and charcoal solutions
- Available with activated charcoal suspension (with or without Sorbitol)
- Twist-Lok™ expressing key aids in charcoal administration
- Easy-to-fill top port
- Available with Universal Solution Spike (for use with pre-filled solution containers)
- No clean up or disinfecting procedures necessary
- All kits come standard with Bite Block

## WITHOUT CHARCOAL

ITEM#	DESCRIPTION	PK
2054	Tum-E-Vac + 3L Solution Bag + 32Fr Lavage Tubing	10
2055	Tum-E-Vac + 3L Solution Bag	10
2075	Tum-E-Vac + Universal Solution Spike	10

## WITH ACTIVATED CHARCOAL

ITEM#	DESCRIPTION	PK
2057	Tum-E-Vac + 32Fr Lavage Tubing + (2) 25g Activated Charcoal	10
2058	Tum-E-Vac + 32Fr Lavage Tubing + (2) 25g Activated Charcoal with Sorbitol	10
2059	Tum-E-Vac + 32Fr Lavage Tubing + (1) 25g Activated Charcoal with Sorbitol + (1) 25g Activated Charcoal	10
2081	Tum-E-Vac + Universal Solution Spike + (2) 25g Activated Charcoal with Sorbitol	10