

PROTECTIVE LAMINATED ULTRASONICALLY WELDED GOWNS

LEVEL 3
NON-STERILE
ULTRASONICALLY WELDED

Our **Level 3 Disposable Protective Laminated Gowns** are impermeable to fluids, offering outstanding protection.

- Fabric Compliant to **EN 20811** and **EN ISO 22610**
- Tested and compliant to **AAMI PB70 - Level 3 Liquid Barrier Performance**
- Tested and compliant to **ASTM F739-12**
- 56 gr/m² soft non-woven spunbond fabric coated with a polyethylene film provides maximum **protection, breathability, and comfort**
- Knitted soft wrist cuffs
- Tie closure
- Ultrasonically welded seams
- Non-Sterile



PRODUCT SELECTION

*ITEM #	**ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
SF-502-S	SF-602-S	S	110 cm / 43"	55 cm / 22"	140 cm / 55"
SF-502-M	SF-602-M	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
SF-502-L	SF-602-L	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
SF-502-XL	SF-602-XL	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
SF-502-XXL	SF-602-XXL	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"
SF-502-XXXL	SF-602-XXXL	XXXL	145 cm / 57"	66 cm / 26"	160 cm / 63"

* PACKAGING SF-502: 40 INDIVIDUALLY BAGGED GOWNS / BOX

** PACKAGING SF-602: 50 BULK PACKED GOWNS / BOX