

# Warm Gel®

## Infant Heel Warmer

Uniform heat to reduce heel sticks

WarmGel® Infant Heel Warmer is a gel-filled, disposable heel warmer designed to provide gentle heat to increase blood flow and enhance infant blood sample collection.

## WarmGel® Technology

### Effective

- Conforms to infant's shape, keeping warmth where needed
- Even heat distribution increases blood flow

### Safe

- Non-toxic food-grade gel won't harm infant's skin or eyes
- Single-patient use minimizes risk of infection

### Reliable

- Reaches maximum temperature in less than 60 seconds, maintains optimal 104° F for 10 minutes <sup>1</sup>
- Trigger disc minimizes auto-activation during shipment

## WarmGel® Versus Liquid

- Conforms to infant's shape
- Provides even heat distribution
- Resists migration, keeping warmth where needed
- Competitively priced compared to liquid-filled products



| Item # | Description                           |
|--------|---------------------------------------|
| 24401  | Liquid Heel Warmer (100/Box)          |
| 20418  | WarmGel® Infant Heel Warmer (100/Box) |

<sup>1</sup>"Neonatal Thermoregulation: Guidelines for Practice," National Association of Neonatal Nurses, 1997.