

Warm Gel®

Get relief where you need it when you need it



Features and Benefits

- Pliable WarmGel® conforms to any body part
- Easy to activate
- Soft non-woven fabric for patient comfort
- Consistent temperature
- Available in two convenient sizes

WarmGel® is a versatile and reliable source for instant heat

Many applications in:

- Maternity
- IV Therapy
- Labor & Delivery
- Rehabilitation

Temperature Range:

WarmGel® has been designed to achieve 112°F (44°C) when activated at room temperature of 75 °F (24 °C)

Item #	Description
20419	Medium 6" x 6" Adult WarmGel® Instant Heat Pack - 40/Case
20420	Large 7" x 9" Adult WarmGel® Instant Heat Pack - 24/Case

