



BATRIK Suture Clamp Tags are designed to protect and safely hold sutures without damage.

- packed in a convenient foam base with easy peel-off (even with gloves) adhesive backing
- the foam base is affixed to a sterile area with the clamps presented in front of the Clamp Tags to allow for easy insertion
- available **sterile** and **non sterile**
- available in 2, 3 or 5 pairs
- radio opaque
- non-toxic
- pre-counted
- latex free



Easy Peel Corner



Suture Clamp Tags

Sterile item #	Non Sterile item #	Packaging	Size	Color
11-1002	911-1002	2 pairs / foam base	Standard \emptyset 0.06" x 0.5" L	Yellow
11-1003	911-1003	3 pairs / foam base	Standard \emptyset 0.06" x 0.5" L	Yellow
11-1005	911-1005	5 pairs / foam base	Standard \emptyset 0.06" x 0.5" L	Yellow
11-1025	911-1025	5 pairs / foam base	Standard \emptyset 0.06" x 0.5" L	Assorted
11-1045	911-1045	5 pairs / foam base	Standard \emptyset 0.06" x 0.5" L	White
11-1065	911-1065	5 pairs / foam base	Standard \emptyset 0.06" x 0.5" L	Blue
11-1085	911-1085	5 pairs / foam base	Standard \emptyset 0.06" x 0.5" L	Red
11-3002	911-3002	2 pairs / foam base	Large \emptyset 0.09" x 0.75" L	Pink
11-3003	911-3003	3 pairs / foam base	Large \emptyset 0.09" x 0.75" L	Pink
11-3005	911-3005	5 pairs / foam base	Large \emptyset 0.09" x 0.75" L	Pink

Our "Easy Peel Corner" is a feature unique to **BATRIK**

Packaging: Sterile: 1 foam base/pouch – 5 pouches/box
Non Sterile: sold per foam base