


BATRIK Clamp Jaw Covers reduce trauma to vessels requiring occlusion, particularly when a hemostatic clamp is utilized. The thread-like indents allow for gentle, constant and non-slip tension.

- manufactured from Medical-Grade silicone
- gamma sterilized
- available in 3 colors and 3 widths to fit most clamp sizes
- convert almost any instrument to atraumatic
- radio opaque
- latex free



Clamp Jaw Covers

Sterile Item #	Non Sterile Item #	Dimensions – inches / cm	Color
12-1000	912-1000	0.1" x 5" / 0.25 x 12.7cm	Red
12-1001	912-1001	0.2" x 5" / 0.51 x 12.7cm	White
12-1002	912-1002	0.3" x 5" / 0.76 x 12.7cm	Blue

Packaging: Sterile: 2 clamps per sterile pouch – 10 pouches per box
Non Sterile: 20 clamps per bag (also available 100 units / bag)