

BATRIK Vascular Silicone Ties permit safe and smooth retraction, occlusion and identification of vessels during surgical procedures.

Unlike clamps and retractors which can pinch the vessels flat, or even crack them, the BATRIK elliptical shape allows for the gentle handling of arteries, veins, nerves, tendons, spermatic cords, ureters, or whenever delicate retraction / occlusion is required.

- manufactured of Medical-Grade silicone
- do not rotate or twist on themselves when unpacked
- easy to separate — eliminating the need for cardboard inserts
- radio opaque
- available in 4 colors and 4 sizes
- latex free



Vascular Silicone Ties

Item #	Dimensions – inches / cm	Color
10-1000	0.075" x 16" / 0.2cm x 40cm	Mini-Red
10-1001	0.075" x 16" / 0.2cm x 40cm	Mini-Blue
10-1002	0.075" x 16" / 0.2cm x 40cm	Mini-Yellow
10-1003	0.075" x 16" / 0.2cm x 40cm	Mini-White
10-1010	0.120" x 16" / 0.3cm x 40cm	Maxi-Red
10-1011	0.120" x 16" / 0.3cm x 40cm	Maxi-Blue
10-1012	0.120" x 16" / 0.3cm x 40cm	Maxi-Yellow
10-1013	0.120" x 16" / 0.3cm x 40cm	Maxi-White
10-1101	0.040" x 16" / 0.1cm x 40cm	Extra Mini-Blue
10-1111	0.200" x 22" / 0.5cm x 56cm	Extra Maxi-Blue
10-1113	0.200" x 22" / 0.5cm x 56cm	Extra Maxi-White

Packaging: 2 ties per sterile pouch — 10 pouches per box