

BATRIK Pin Guards You will no longer have to use tape or tubing as our BATRIK Pin Guards will safely cover the sharp ends of your K-Wire and Steinmann Pins.

- color coded by wire size which eliminates the selection of the wrong size
- a removable sleeve is provided to avoid injury and facilitate application (sterile only)
- packaging designed to reduce waste
- non sterile (100 units / bag) available with or without plastic sleeves
- latex free



Pin Guards

Sterile item #	Non Sterile item # ⁽¹⁾	Color	Dimensions inches / mm	Sterile Packaging	Non-Sterile Packaging**
14-1028	914-1028	Yellow	0.028" / 0.71mm	80 units/box*	100 units/box
14-1035	914-1035	Light Blue	0.035" / 0.89mm	80 units/box*	100 units/box
14-1045	914-1045	White	0.045" / 1.14mm	80 units/box*	100 units/box
14-1054	914-1054	Dark Red	0.054" / 1.40mm	80 units/box*	100 units/box
14-1062	914-1062	Green	0.062" / 1.57mm	80 units/box*	100 units/box
14-1078	914-1078	Pink	0.078" / 1.98mm	80 units/box*	100 units/box
14-1093	914-1093	Grey	0.093" / 2.36mm	80 units/box*	100 units/box
14-1109	914-1109	Beige	0.109" / 2.77mm	80 units/box*	100 units/box
14-1125	914-1125	Purple	0.125" / 3.17mm	80 units/box*	100 units/box
14-1028-1	N/A	Yellow	0.028" / 0.71mm	40 units/box	N/A
14-1035-1	N/A	Light Blue	0.035" / 0.89mm	40 units/box	N/A
14-1045-1	N/A	White	0.045" / 1.14mm	40 units/box	N/A
14-1062-1	N/A	Green	0.062" / 1.57mm	40 units/box	N/A
14-1141	914-1141	Brown	0.141" / 3.58mm	80 units/box*	100 units/box
14-1156	914-1156	Dark Blue	0.156" / 3.96mm	80 units/box*	100 units/box
14-1188	914-1188	Red	0.188" / 4.77mm	80 units/box*	100 units/box

*Sterile Packaging: 2 units per sterile pouch; 40 pouches per box.

** Non-Sterile: Does not include plastic sleeve

⁽¹⁾ Non-sterile sold without plastic sleeves. See price list for part number and price with sleeve.