


BATRIK Magnetic Mats provide a **sterile, neutral field** for the safe placement and **hands-free transfer** of instruments.

- **raised magnets** hold instruments firmly in place while allowing easy grasping of instruments
- **round** magnets to avoid silicone rubber perforations
- **flexible** and **non-slip** backing drape conforms to any surface shape

Reusable and Autoclavable (up to 134°C)

Latex free



Magnetic Mats

Item #	Dimensions – inches / cm	Packaging
MAG-02-1	12" x 15" / 30cm x 38cm	1 unit
MAG-02-5	12" x 15" / 30cm x 38cm	5 units