


### BATRIK Endoscopy BO2™ Bite Blocks with Nasal Oxygen Delivery

significantly enhances patient safety by providing supplemental oxygenation during gastrointestinal endoscopy preventing hypoxia in sedated patients, while ensuring superior patient comfort.

- lower procedural costs by obviating the need for separate nasal prongs and oxygen masks, the BO2 Bite Block helps save time and money
- ensure patient comfort and safety as with other Bite Blocks from BATRIK Medical, the BO2 Bite Block is designed to be free of sharp edges to ensure patient comfort and safety. Additionally, our bite blocks are guaranteed to be latex-free, lead-free and DEHP-free.
- choose from disposable and reusable options. The BO2 Bite Block is available in single use and reusable options. Use our flexible ARC Wire Brush #41-32505, which will allow you to properly clean the nasal prongs and inlet of the reusable BO2.
- ensure secure fit. An optional retention strap is available to ensure a secure fit. These straps, like the bite block, are latex free, lead free and DEHP free.




#### BO2™ Bite Block

Item #	Description	Packaging
BO2	18mm / 54FR - Single use, includes strap	100 / Box
BO2-N	18mm / 54FR - Single use, without strap	100 / Box
BO2-RN	18mm / 54FR - Reusable, without strap	25 / Box
BO2-R	18mm / 54FR - Reusable, includes strap	25 / Box

### BATRIK Endoscopy Bitestrong™ Bite Blocks

are designed to provide a safe passage through the patient's mouth during TEE thus protecting the patient's teeth as well as the endoscope. Strict quality control procedures and the careful selection of raw materials used in the manufacture of our Bite Blocks is your guarantee that they are latex free, lead free and DEHP free.

#### Single-Use



#### Bite Blocks

Item #	Description	Packaging
B4	19mm / 60FR - Single use	100 / Box
B4T	19mm / 60FR - Single use, includes strap	100 / Box
B4T-N	19mm / 60FR - Single use, without strap	100 / Box
B1-25	16mm / 48FR - <b>Pediatric</b> single use, includes strap	25 / Box
B1	16mm / 48FR - <b>Pediatric</b> single use, includes strap	100 / Box
B1-N	16mm / 48FR - <b>Pediatric</b> single use, without strap	100 / Box

**BEST SELLER**

**BATRIK Endoscopy Bitestrong™ Single Use Bite Blocks**

- our molding process ensures a smooth texture finish with no sharp edges (flashing) that can often irritate a patient’s mouth
- our bite blocks are individually wrapped in sealed, easy-to-open PE bags
- available without dental retention
- additional straps are also available, if required




**Bitestrong™ Single Use Adult Bite Block**

Item #	Description	Packaging
B3	20mm / 60FR - Single use, includes strap	100 / Box
B3-N	20mm / 60FR Reusable, without strap	100 / Box
<b>NEW</b> B3-SA	Same as B3, with strap <u>attached</u> to one side	100 / Box
B3-NDR	Same as B3, no dental retention rim	100 / Box
B3-NDR-N	Same as B3-N, no dental retention rim, without strap	100 / Box

**BATRIK Endoscopy Bitestrong™ Reusable Bite Blocks**

- our reusable Bite Blocks are designed for repeated usage and will withstand commonly practiced chemical and steam sterilization cycles (134°C / 275°F)<sup>(1)</sup>
- our molding process ensures a smooth textured finish with no sharp edges (flashing) that can often irritate a patient’s mouth
- our Bite Blocks are individually wrapped in sealed, easy-to-open PE bags
- available with or without an adjustable retention strap
- additional straps are also available, if required



**Bitestrong™ Reusable Adult Bite Block and Straps**

Item #	Description	Packaging
B5	20mm / 60FR Reusable, includes strap	100 / Box
B5-N	20mm / 60FR Reusable, without strap	100 / Box
B5-25	20mm / 60FR Reusable, includes strap	25 / Box
B5-N-25	20mm / 60FR Reusable, without strap	25 / Box
B6	19mm / 60FR Reusable	100 / Box
B6T	19mm / 60FR Reusable, includes strap	100 / Box
B6T-N	19mm / 60FR Reusable, without strap	100 / Box
B6-25	19mm / 60FR Reusable	25 / Box
B6T-25	19mm / 60FR Reusable, includes strap	25 / Box
B6T-N-25	19mm / 60FR Reusable, without strap	25 / Box
BS	Single-Use Bite-Block straps 55cm x 2cm / 21.7" x 0.79"	100 / Box

Order your single use Bite Block sample box today! **Item # B-SAMPLE-S**  
The BATRIK BITESTRONG BiteBlock Single Use Sample Box includes:

Order your reusable Bite Block sample box today! **Item # B-SAMPLE-R**  
The BATRIK BITESTRONG BiteBlock Reusable Sample Box includes: