

BATRIK Safe-Soak Endo Sponge™ are specifically fabricated for use with all endoscopes. Unlike other colored sponges, our sponges are near white as they will better evidence the effectiveness of the cleaning process. Dark colors easily hide the contaminants, where with a white sponge you can visually evaluate the pre-clean efficacy.

Available in two formats, Flat or Tubular. Both provide the identical advantages:

- 100% Medical Grade High Density Foam
- will not flake or scratch expensive equipment
- high absorbency as it retains more solution than conventional sponges or paper wipes
- compatible with Enzymatic or Disinfectant solutions
- green friendly



Or use it as a bumper to protect your scopes!

Safe-Soak Endo Sponge™

Item #	Description
17-100	Tubular Sponge / White Dispenser Box
17-101	Tubular / Bulk
17-200	Flat / White Dispenser Box
17-201	Flat / Bulk

Packaging
25 Units / Box – 4 Boxes / Case
560 Units / Case
50 Units / Box – 4 Boxes / Case
1000 Units / Case

Note: The white sponge may at times appear slightly yellowed. This is only due to UV exposure, not a reflection to its age.