

BATRIK ENDOSCOPY Single Use ARC™ Channel Cleaning Brushes

BATRIK offers one of the largest selection of Endoscopy Channel Cleaning Brushes in the industry. Our brushes are available in a variety of lengths, material, and styles to accommodate the widest range of endoscopic devices used in the healthcare industry today.

No matter what your requirements, you'll be able to find the brush that suits your specific needs. If you don't see the brush you require, just call us. We have well over 50 different channel brushes available and if you can't find the exact one that you need, no problem — we can custom fabricate to your specifications.

Why Choose a BATRIK ARC Channel Brush?

- one of the industry's longest brush heads on the market: 15mm to 24 mm in length. More bristles translate into more cleaning contact within the channels
- bristle fan tip on the head of the brush protects and cleans the channel. No sharp end or hard metal tips that can scratch or cause damage
- latex free
- available in two different materials in order to accommodate different customer preferences:
 - POLYPROPYLENE is the industry standard as it is a low cost material and offers greater flexibility.
 - DELRIN is slightly more expensive, but due to its more rigid feeling, it is easier to insert and manipulate. This in turn reduces time and labor which can offset the added cost. ALL BATRIK small diameter brushes are supplied only with DELRIN, otherwise they would be too fragile



F= receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.

K= receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

FK= includes both sponge and valve cleaning brush.



5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.85 MM – PUNCHY PINK POLYPROPYLENE

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
311-23050	90" / 230cm	0.94" / 24mm	0.20" / 5mm	0.073" / 1.85mm	Single
311-23050-2	90" / 230cm	0.94" / 24mm	0.20" / 5mm	0.073" / 1.85mm	Double
311-18050	70" / 180cm	0.94" / 24mm	0.20" / 5mm	0.073" / 1.85mm	Single
311-10050	39" / 100cm	0.94" / 24mm	0.20" / 5mm	0.073" / 1.85mm	Single

Packaging: 50 / Box

Please note that you may receive the brush in a RED color instead of PINK. It is perfectly normal as we are slowly phasing out the RED color. The quality of the product remains identical (31 Series).

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.85 MM – PUNCHY PINK POLYPROPYLENE

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
333-23056	90" / 230cm	0.79" / 20mm	0.20, 0.24" / 5, 6mm	0.073" / 1.85mm	Single
333-23056-2	90" / 230cm	0.79" / 20mm	0.20, 0.24" / 5, 6mm	0.073" / 1.85mm	Double

Packaging: 50 / Box

Please note that you may receive the brush in a RED color instead of PINK. It is perfectly normal as we are slowly phasing out the RED color. The quality of the product remains identical (33 Series).

BATRIK ENDOSCOPY Single Use ARC™ Channel Cleaning Brushes**5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM – METALLIC GREY DELRIN**

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
300-18050	70" / 180cm	0.79" / 20mm	0.20" / 5mm	0.070" / 1.75mm	Single
300-23050	90" / 230cm	0.79" / 20mm	0.20" / 5mm	0.070" / 1.75mm	Single
300-23050-2	90" / 230cm	0.79" / 20mm	0.20" / 5mm	0.070" / 1.75mm	Double

Packaging: 50 / Box

Please note that you may receive the brush in a **WHITE** color instead of **GREY**. It is perfectly normal as we are slowly phasing out the **WHITE** color. The quality of the product remains identical (30 Series)

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.75 MM – METALLIC GREY DELRIN

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
322-23056	90" / 230cm	0.79" / 20mm	0.20, 0.24" / 5, 6mm	0.070" / 1.75mm	Single
322-23056-2	90" / 230cm	0.79" / 20mm	0.20, 0.24" / 5, 6mm	0.070" / 1.75mm	Double

Packaging: 50 / Box

Please note that you may receive the brush in a **WHITE** color instead of **GREY**. It is perfectly normal as we are slowly phasing out the **WHITE** color. The quality of the product remains identical (32 Series)

**3 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.55 MM – METALLIC PURPLE DELRIN**

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
344-23030	90" / 230cm	0.79" / 20mm	0.12" / 3mm	0.060" / 1.55mm	Single
344-18030	70" / 180cm	0.79" / 20mm	0.12" / 3mm	0.060" / 1.55mm	Single
344-10030	39" / 100cm	0.79" / 20mm	0.12" / 3mm	0.060" / 1.55mm	Single

Packaging: 50 / Box

Please note that you may receive the brush in a **GREEN** color instead of **PURPLE**. It is perfectly normal as we are slowly phasing out the **GREEN** color. The quality of the product remains identical (34 Series)

**2.5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM – YELLOW DELRIN**

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
35-23025	90" / 230cm	0.59" / 15mm	0.098" / 2.5mm	0.039" / 1.00mm	Single

Packaging: 50 / Box**2 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM – YELLOW DELRIN**

Item #	Length inches / cm	Brush Length inches / mm	Brush (Ø) Inches / mm	Sheath (Ø) Inches / mm	Brush Ends
35-23020	90" / 230cm	0.59" / 15mm	0.079" / 2.0mm	0.039" / 1.00mm	Single

Packaging: 50 / Box