

BATRIK Valve / Control Cleaning Brushes

VALVE / CONTROL CLEANING BRUSH – 5/7



Item #	Length inches / cm	Brush Length inches / mm	Brush 1 (Ø) inches / mm	Brush 2 (Ø) inches / mm	Description
27-05197	5.00" / 13 cm	1.00" / 25 mm	0.20" / 5 mm	0.28" / 7 mm	Dble. Trumpet

Packaging: (Single-use recommended) Bag of 3 - **ITEM #** 27-05197-50 / Bag of 50

VALVE / CONTROL CLEANING BRUSH – 5/11



Item #	Length inches / cm	Brush 1 Length inches / mm	Brush 2 Length inches / mm	Brush 1 (Ø) inches / mm	Brush 2 (Ø) inches / mm	Description
27-57787	6.00" / 15cm	1.18" / 30 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm	Dble. Universal

Packaging: (Single-use recommended) Bag of 3 - **ITEM #** 27-57787-50 / Bag of 50

VALVE / CONTROL CLEANING REUSABLE BRUSH – 5/11 – EXTRA LONG



Item #	Length inches / cm	Brush 1 Length inches / mm	Brush 2 Length inches / mm	Brush 1 (Ø) inches / mm	Brush 2 (Ø) inches / mm	Description
40-18511	8.07" / 20.5cm	1.57" / 40 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm	Dble. Universal

Packaging: (Re-usable and can be steam sterilized) Bag of 3

VALVE / CONTROL CLEANING BRUSH – SINGLE USE BRUSH



or



Item #	Length Inches / cm	Brush Length Inches / mm	Brush (Ø) Inches / mm	Description
36-10010*	4.3" / 11 cm	1.57" / 40 mm	0.39" / 10 mm	Single end large

Packaging: (Single-use recommended) Bag of 3 - **ITEM #** 36-10010-50 / Bag of 50

*Actual brush handle may vary, but brush head is unchanged.