

BATRIK ENDOSCOPY V-Brush - Reusable External Cleaning Brush External cleaning brush for endoscopes / optical instruments. With its unique shape and innovative design, the BATRIK Reusable External Cleaning Brush is the ultimate tool for cleaning the exterior of your endoscopic instruments without the risk of damaging expensive equipment.

- easy-grip handle ensures comfort and control
- numerous dense bristles provide additional surface contact
- superior endoscopic scope cleaning with a brush specifically designed for this purpose



Handle Length / Diameter
120mm / 12mm

Brush Length / Diameter
90mm / 25mm

Total Length
210mm

- polyamide (TYNEX®) bristles
- superior bristle durability
- steamable thermoplastic handle
- ergonomic handle
- steam / chemical sterilizable
- latex free

Reusable External Cleaning Brush

Item #	Description	Packaging
36-20010	Reusable External Cleaning Brush	3 units / Bag

BATRIK ENDOSCOPY V-Brush Single Use External Cleaning Brush External cleaning brush for endoscopes / optical instruments. An ideal fit for most endoscopic instruments, the BATRIK Single Use External Cleaning Brush will provide excellent cleaning capabilities. It's U-shape bristle design provides continuous and effective cleaning. There are no sharp or metallic ends that can scratch or cause damage to valuable equipment.

Many of the same features as the reusable brush

- polypropylene bristles
- vinyl caps
- ergonomic handle
- latex free



Handle Length / Diameter
80mm / 3mm

Brush Length / Diameter
70mm / 25mm

Total Length
150mm

Single Use External Cleaning Brush

Item #	Description	Packaging
36-20005	Single Use External Cleaning Brush	50 units / Box