


BATRIK Corner Pockets™ Avoid frustration and costs from having to reprocess your surgical trays due to ripped sterile wrapping. Corner Pockets™ will prevent tearing associated with sharp corners or feet on surgical trays. They will fit all types of corners (even the round ones), and the bright orange-red color makes them easily identifiable and will reduce the risk of being discarded. We are also pleased to introduce one of the only non-foam single-use corner protectors (see below).

- no moisture retention
- reusable or single-use options
- non sterile
- latex free



The Reusable corner pockets can withstand repeated sterilization up to 134°C; after 10 cycles, no product degradation was observed.

Corner Pockets™

| Item # | Description | Packaging |
|-----------|------------------------------|------------------------|
| 18-CP07 | Corner Pockets™ — Reusable | 80 pcs / Bag |
| 18-CORNER | Corner Pockets™ — Reusable | 4 Bags / Box (320 pcs) |
| 18-EVAL7 | Corner Pockets™ — Reusable | 24 pcs / Bag |
| 18-CP02* | Corner Pockets™ — Single-Use | 240 pcs / Bag |
| 18-EVAL2* | Corner Pockets™ — Single-Use | 24 pcs / Bag |



AN INDUSTRY FIRST – SINGLE USE CORNER POCKETS!

Our Single-Use corner pockets have the following advantages:

- save time; **easy** to use and position
- **no moisture** retention


** Although the product appears similar to the reusable device, please be advised that it is blended differently to reduce costs and has only been validated as single use. Re-using the single use Corner Pocket may compromise sterility.*

***NEW* FOAM CORNER PROTECTORS**

Soft Pockets™

| Item # | Description | Packaging |
|-----------|---|---------------|
| 18-CPF-02 | Soft Pockets™ - Single-Use Foam Corner Protectors | 80 sets / Box |

“Our foam is 0.25” thicker than the competitors. Especially useful with heavier Instrument sets. Our size difference is hardly noticeable, but the benefits are very apparent.”
(24 boxes per pallet)



6” x 6” x 1.50”
(15.2cm x 15.2cm x 3.8cm)