

BATRIK Instrument Protectors will protect instrument tips and edges from damage during handling and sterilization.


- available in a wide array of sizes for almost every application
- sizes are color coded
- opaque and tinted protectors are radiopaque
- steam and gas permeable
- non sterile
- latex free

Variations Include

- perforated protectors (holes) which reduce vacuum seal and moisture build-up
- osteotome protectors
- double protectors
- retractor protectors
 - “Weitlaner” Perforated “Square” Style
 - Single-Tip Perforated Clear “Tube” Style
- validated for pre-vacuum sterilization



Tinted



Opaque – Solid Color

Item # (regular)	Item # (perforated)
1 RS-01	PS-01
2 RS-02	PS-02
3 RS-03	PS-03
4 RS-04	PS-04
5 RS-05	PS-05
6 RS-06	PS-06
7 RS-07	PS-07
8 RS-08	PS-08
9 RS-09	PS-09
RS-10	PS-10

Packaging: 100 units/Bag

Tinted Color

Item # (regular)	Item # (perforated)	Color
RT-01	PT-01	White
RT-02	PT-02	Blue
RT-03	PT-03	Green
RT-04	PT-04	Red
RT-05	PT-05	Yellow
RT-06	PT-06	Maroon
RT-07	PT-07	Orange
RT-08	PT-08	Brown
RT-09	PT-09	Black
RT-10	PT-10	Assorted

Packaging: 100 units/Bag

Clear - No Color

Item # (regular)	Item # (perforated)
1 RC-01	PC-01
2 RC-02	PC-02
3 RC-03	PC-03
4 RC-04	PC-04
5 RC-05	PC-05
6 RC-06	PC-06
7 RC-07	PC-07
8 RC-08	PC-08
9 RC-09	PC-09
RC-10	PC-10

Packaging: 100 units/Bag



Clear


Solid

Double Protectors

Item #	Dimensions – inches / mm	Color
DG-01	0.0625 x 0.75" / 1.59 x 19.05mm	White
DG-02	0.078 x 0.75" / 1.98 x 19.05mm	Blue
DG-03	0.109 x 0.75" / 2.78 x 19.05mm	Green
DG-04	0.125 x 1" / 3.18 x 25.4mm	Red
DG-05	0.1875 x 1" / 4.76 x 25.4mm	Yellow
DG-07	0.0625 x 0.375 x 1" / 1.59 x 9.53 x 25.4mm	Orange
DG-10	Assorted	Assorted



Packaging: 50 units/Bag - Perforated only



Osteotome Protectors

Item #	Dimensions – inches / mm
OP-11	0.25 x 0.75" / 6.35 x 19.05mm
OP-12	0.375 x 1" / 9.53 x 25.4mm
OP-13	0.5 x 1" / 12.7 x 25.4mm
OP-14	0.75 x 1" / 19.05 x 25.4mm
OP-15	1 x 1" / 25.4 x 25.4mm
OP-16	1.25 x 1.25" / 31.75 x 31.75mm
OP-17	1.5 x 1.5" / 38.1 x 38.1mm
OP-18	1.75 x 1.25" / 44.45 x 31.75mm
OP-20	Assorted

Packaging: 50 units/Bag - Perforated only



Retractor Protectors "Weitlaner"

Item #	Size	Color
RGM-01	Small	Pink
RGM-02	Small	Clear
RGM-03	Medium	Green
RGM-04	Medium	Clear

Packaging: 50 units/Bag

Retractor Protectors Single-Tip

Item #	Dimensions – inches / cm	Packaging
RGS-01	1.5" / 3.8cm	55 units / Bag
RGS-02	2" / 5.1cm	45 units / Bag
RGS-03	3" / 7.62cm	40 units / Bag